



# REVISION

## KNEE SYSTEM

ADAPTIVE • SMART • COMPLETE

SURGICAL TECHNIQUE  
SURGICAL TECHNIQUE



# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Index

|                                             |                |
|---------------------------------------------|----------------|
| Indications, Contraindications and Warnings | Pag. >> 10     |
| Introduction                                | Pag. >> 13     |
| <br>SURGICAL TECHNIQUE                      | <br>Pag. >> 14 |
| CCK Tibia                                   | Pag. >> 15     |
| CCK Femur                                   | Pag. >> 24     |
| CCK Trial reduction                         | Pag. >> 33     |
| CCK Definitive implant                      | Pag. >> 35     |
| H Tibia                                     | Pag. >> 41     |
| H Femur                                     | Pag. >> 49     |
| H Trial reduction                           | Pag. >> 58     |
| H Definitive implant                        | Pag. >> 60     |
| Extraction of Hinge components              | Pag. >> 69     |
| Patellar prosthesis                         | Pag. >> 71     |
| <br>INSTRUMENT SET                          | <br>Pag. >> 73 |
| Instruments needed for Hinge only           | Pag. >> 82     |
| <br>PRODUCT COMBINATIONS                    | <br>Pag. >> 83 |
| <br>PRODUCT CODES                           | <br>Pag. >> 84 |
| <br>TABLES SUMMARY                          | <br>Pag. >> 92 |

*Limacorporate S.p.A. is a manufacturer of prosthetic implants and as such does not perform medical procedures.*

*This documentation concerning surgical techniques, which provides surgeons with general guidelines for implanting the REVISION KNEE SYSTEM, was developed with the advice of a team of surgical experts. All decisions as to the type of surgery and most suitable technique are obviously the responsibility of the health care professional. Surgeons must make their own decisions as to the adequacy of each planned implant technique based on their training, experience and the clinical condition of the patient.*


*For further information about our products, please visit our web site at [www.limacorporate.com](http://www.limacorporate.com)*





# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Surgical Steps






Articular test



CCK or H Definitive femoral implant assembly



CCK definitive Tibial implant assembly




CCK definitive Tibial insert assembly



**MULTIGEN** CCK definitive Implant



H definitive Tibial implant assembly



H Hinge mechanism final assembly



**MULTIGEN** H definitive Implant




# REVISION

KNEE SYSTEM

ADAPTIVE • SMART • COMPLETE





## THE MULTIGEN-PLUS H TOGETHER WITH THE MULTIGEN-PLUS CCK, CONSTITUTE A COMPLETE REVISION KNEE SYSTEM

### REGAINING MOBILITY

The LimaCorporate Multigen Plus Total Knee system provides the surgeon with a solution to most knee surgery cases, from simple to complex primary and revision.

## MODULAR

Different modules to adapt the stem positioning and improve the fit for each individual patient.

## SMART

The interchangeability of the components and a common instrument set for the CCK and the Hinge facilitates the surgical steps.

## COMPLETE

The system can solve several challenging situations in revision surgery thanks to a condylar constrained and a hinged knee configuration, designed to maximize stability alongside mobility.

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Indications, Contraindications and Warnings

### ▼ INDICATIONS



Please follow the instructions for use enclosed in the product packaging.

- Advanced articular destruction generated by primary degenerative or post-traumatic arthrosis or rheumatoid arthritis;
- Traumatic events with articular pain;
- Avascular necrosis;
- congenital or acquired deformity;
- failures of previous operations, articular reconstruction, arthrodesis, hemi- arthroplasty or total arthroplasty.

### ▼ CONTRAINDICATIONS

- Acute or chronic infections, local or systemic infections;
- serious muscular, neurological or vascular diseases affecting the concerned limb;
- bone destruction or poor bone quality, which could compromise implant's stability;
- any concomitant disease and dependence that might affect the implanted prosthesis;
- allergy to material.

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Indications, Contraindications and Warnings

### ▼ RISK FACTORS

The following risk factors may lead to poor results by implanting this prosthesis:

- overweight - According to the definition of the World Health Organization (WHO), Body Mass Index (BMI) greater than or equal to 25 Kg/m<sup>2</sup>;
- strenuous physical activities (active sports, heavy physical work);
- fretting of modular junctions;
- incorrect implant positioning;
- insufficient bone to support the femoral and/or tibial components;
- medical disabilities which can lead to an unnatural gait and loading of the knee joint;
- muscle deficiencies;
- multiple joint disabilities;
- refusal to modify postoperative physical activities;
- patient's history of infections or falls;
- systemic diseases and metabolic disorders;
- local or disseminated neoplastic diseases;
- drug therapies that adversely affect bone quality, healing, or resistance to infection;
- drug use or alcoholism;
- marked osteoporosis or osteomalacia;
- patient's resistance to disease generally weakened (HIV, tumour, infections);
- severe deformity leading to impaired anchorage or improper positioning of implants;
- mistakes while performing the surgical technique;
- Lack of tightening of the hinge locking screw (only for MULTIGEN PLUS H).

**NOTE.** With CCK-H it is always necessary to use modules in combinations with stems.

### ▼ WARNINGS

Surgeons must carefully plan the procedure after viewing the surgical technique to be used to implant this system.

#### ALLOWED/PROHIBITED COMBINATIONS

For the MULTIGEN-PLUS CCK system, the following combinations of sizes are allowed:

| MULTIGEN <sup>+</sup> CCK   |    | FEMORAL COMPONENT |    |    |    |    |
|-----------------------------|----|-------------------|----|----|----|----|
|                             |    | #1                | #2 | #3 | #4 | #5 |
| TIBIAL COMPONENT= TIB. LIN. | #1 | OK                | OK | NO | NO | NO |
|                             | #2 | OK                | OK | OK | NO | NO |
|                             | #3 | NO                | OK | OK | OK | NO |
|                             | #4 | NO                | NO | OK | OK | OK |
|                             | #5 | NO                | NO | NO | OK | OK |

For the MULTIGEN PLUS H system, the following combinations of sizes are allowed:

| MULTIGEN <sup>+</sup> H |    | FEM. COMP.= TIB. LIN. |    |    |    |    |
|-------------------------|----|-----------------------|----|----|----|----|
|                         |    | #1                    | #2 | #3 | #4 | #5 |
| TIBIAL COMPONENT        | #1 | OK                    | OK | NO | NO | NO |
|                         | #2 | OK                    | OK | OK | NO | NO |
|                         | #3 | OK                    | OK | OK | OK | NO |
|                         | #4 | OK                    | OK | OK | OK | OK |
|                         | #5 | OK                    | OK | OK | OK | OK |

- The CCK Tibial modules can be assembled only with CCK tibial component and with a CCK/H stem.
- The Tibial augments can be assembled only with cemented fixed tibial component.
- The H Tibial modules can be assembled only with an H tibial component and with a CCK/H stem.
- The H Tibial augments can be assembled only with an H tibial component.
- The H Tibial liners must be coupled only with H tibial component.



### ▼ INTRODUCTION


In Revision Knee Surgery (Revision TKA) the surgeon removes a previous implant and replaces it with a new prosthesis. Revision implants may also be used in primary cases with severe bone loss, traumatic events or compromised ligaments.

Failures in primary TKA are not common, but actually they result to be very problematic for the patient. TKA patients in need of a revision, have often compromised collateral ligaments and subsequently suffer from varus/valgus instability. The common device used to treat these cases is a Constrained Condylar Knee (CCK) implant.

Rotating-Hinge total knee prostheses (H) is used for the treatment of global instability due to weakness or absence of collateral ligaments or severe bone loss around the joint. A ligament deficiency in the knee can be related to a previous trauma, associated with severe varus-valgus deformities, or due to revision surgery associated with severe bone losses. The rotating-hinge arthroplasty offers enough stability, allowing an intrinsic rotation that stimulates the biomechanical reply of a normal knee and reduces the stress produced by an elevated constriction.

With the REVISION KNEE SYSTEM, MULTIGEN-PLUS CONDYLAR CONSTRAINED KNEE together with the MULTIGEN-PLUS HINGE KNEE, Lima Corporate, provides a solution for most of knee surgery situations from complex primary to revision cases.

**NOTE.** *The Multigen plus system has the same A/P resections on the femoral component for primary and revision therefore it is a bone sparing system.*



# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Surgical Technique




Figure 1



Figure 2

### ▼ APPROACH

When possible, follow the scar left by the primary procedure (*Fig. 1*). Where parallel incisions are present, the more lateral is usually preferred, as the blood supply to the extensor surface is medially dominant.

### ▼ CAPSULAR INCISION

The fascial incision extends from the rectus femoris proximal margin to the distal margin of the tibial tubercle following the patella's medial border.

Where mobilisation of the extensor mechanism and patella is problematic, extend the skin and capsular incisions proximally.

### NOTE.

*CCK or Hinged knee can be used:*

- in a revision surgery to replace a failed primary system depending on the ligaments;
- in cases of bone defects, trauma, tumors or other non standard situations ( In this case there is no necessity of implant removal).

### ▼ PREVIOUS IMPLANT REMOVAL

Remove the previously implanted components as well as any residual cement from bone surfaces, taking care to preserve as much bone as possible (*Fig. 2*).



Figure 3



Figure 4

## ▼ TIBIA

### TIBIAL STEM SIZE DETERMINATION

Introduce the first smallest reamer (14mm) into the tibial canal (Figs. 3-4).

Progressively ream the canal until cortical contact is achieved.

If the stability is reached at the reference number S, it means that a short stem is needed (Fig. 5).

The trial short stem (9066.47.230 - 9066.47.255) needs to be assembled in combination with the correct tibial offset trial module (0/+3/+6).

For the definitive implant the correct stem diameter + SHORT tibial module straight +/-3mm/+6mm should be taken.

**NOTE.** The thickness of the stem should always be the thinnest one when in between 2 sizes, otherwise can lead to stem tip pain.




Figure 5

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### CCK Tibia




Figure 6

If the second line (L) is reached, it means a long stem is needed (*Fig. 6*). Therefore use the trial tibial stem connected with the chosen tibial trial offset module (straight/+3mm/+6mm) and the extension module (long) (9066.47.615). The final implant will consist of the definitive stem and the “L” module (*Tab. 1*).

Tab. 1

| Ref      | LENGTH      | TRIAL                            | DEFINITIVE      |
|----------|-------------|----------------------------------|-----------------|
| <b>S</b> | <b>30mm</b> | S Tibial Module                  | S Tibial Module |
| <b>L</b> | <b>55mm</b> | S Tibial Module +<br>L Extension | L Tibial Module |

Whether the reamer goes below the L reference use reamers with the highest reference point as diameters.

If there is a difference between the medial and the lateral plateau use the highest as reference.

The size of the latest reamer used corresponds to the diameter of both trial and definitive stem.




Figure 7

### TIBIAL RESECTION (CCK TIBIAL SLOPE 6°)

Remove the handle from the reamer and insert the prismatic guide onto the reamer shaft. Insert the tibial cutting guide (9066.22.120) on the prismatic guide.

**NOTE.** The CCK tibia resection guide has a 6° of tibial slope (9066.22.120) (Fig. 7). Please be aware of the correct rotation.

Rotate the millimetric screw until the "0" marking is reached on the tibial cutting guide.

Introduce the sickle (9066.12.010) or a free blade with thickness 1.27mm through the cutting guide slot and place it on the highest tibial plateau (Fig. 7) (if there is any difference between the medial and the lateral plateau).

Lock the screw of the prismatic guide on the reamer shaft. Rotate the millimetric screw until the desired resection level is reached (Fig. 8) and secure the tibia resection guide using pins (Fig. 9).

It is also possible to use the +0mm tibia depth resection gauge adapting the space (9066.25.160) to measure the gap between the two tibial emiplateau.




Figure 8




Figure 9

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## CCK Tibia



Figure 10

Proceed with the tibial resection (Fig. 10).

Whether a tibial augment is needed it's possible to directly perform the +7mm or the +12mm resection to host it. Rotate the millimetric screw until the desired resection level +7mm or +12mm, referring to the last performed resection, is reached (Fig. 11).

Remove the pins, the IM tibial resection guide and the prismatic guide.

If needed, it is possible to perform the tibial augment cuts after the tibial resection guide has been removed, using the augment tibial resection guide (9066.47.055) that has to be inserted onto the reamer shaft as well (Fig. 12).

You can also use the standard IM tibial cutting guide, where there is a need for augments (Figs. 13-14).

Whether an augment is needed it's possible to directly perform the +7mm or the +12mm resection through the dedicated slot on the IM Tibial cutting guide for CCK (9066.47.050)





Figure 11



9066.47.050

Figure 13




Figure 12



Figure 14

2 different trial tibial augments can be used:

7mm (9066.47.315/325/335/345/355)

12mm (9066.47.320/330/340/350/360)




Figure 15

### TIBIAL PLATE SIZE AND TIBIAL STEM OFFSET DETERMINATION

Leaving the reamer in situ, place the trial tibial plate (9066.47.050-150) on the tibial surface.

Slide the straight tibial dialer (9066.47.170) along the arm and insert it into the circular slot of the trial tibial plate.

At the straight tibial dialer corresponds the trial short straight tibial module (9066.47.200) (without offset).

Select the tibial size to achieve maximal tibial coverage (Figs. 15-16).

**NOTE.** In a CCK surgery be aware of using the tibial trial plates without the marking "H".




Figure 16

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### CCK Tibia



Figure 17




Figure 18

If a tibial augment is required, a trial augment (9066.47.315-360) of the same size of the tibial plate can be fixed magnetically to the trial plate.

To obtain a better coverage of the resected tibial surface, it could be necessary to offset the trial tibia plate replacing the neutral dialer with the +3mm or +6mm eccentric dialer and rotating until the best positioning is reached (Fig. 17).

A +3mm and +6mm dialer corresponds to a +3mm and +6mm tibial module (see Tab. 2).

Tab. 2

| Dialer | Trial/Definitive                         |
|--------|------------------------------------------|
| 0      | trial short straight tibial module       |
| +3     | +3 trial short offset +3mm tibial module |
| +6     | +6 trial short offset +6mm tibial module |

**NOTE.** Once the right positioning is reached, take note of the number marked on the dialer aligned with the reference on the trial tibial plate (Fig. 18).



Figure 19

### TIBIAL KEEL PREPARATION

Insert two K-wires (dia. 2mm) through the small diameter holes of the CCK trial tibial plate (Fig. 19).

**NOTE.** The K-wires are not included in the set.

Remove the dialer, the CCK trial tibial plate and the reamer.

Reposition the CCK trial tibial plate using the two K-wires as guides (Fig. 20) and secure the tibial plate using the pins.

If a trial tibial augment is used, secure the plate using the long pins.

Remove the two K-wires.

**NOTE.** Make sure to use the winged broach for CCK (without H) (9066.47.180).

Insert the guide for the winged broach onto the trial tibial plate (Fig. 21).




Figure 20




Figure 21

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## CCK Tibia

Insert the CCK winged broach into the guide and hammer it axially until the lever tooth engages the tibial winged broach to the guide for the winged broach (Fig. 22).

Use the extractor to remove the broach and the guide (Fig. 23).

In order to avoid any interference between tibial module and bone, reattach the guide for the broach onto the trial tibial plate and ream the tibia using the tibial reamer for offset until it's in even contact with the guide for the tibial broach (Fig. 24).

This procedure is necessary when a tibial module with +6mm offset needs to be used.

If a tibial module with + 3mm offset is used this procedure is necessary only if the last reamer was smaller than 20mm.

To check the tibial alignment, the alignment rod (9066.15.090) can be used through the slot of the handle (Fig. 24).



Figure 22

Figure 23



Figure 24



Figure 25




Figure 26



Figure 27

### TIBIAL TRIALS

Connect the CCK trial tibial plate with the CCK tibial keel (9066.47.180) and the chosen tibial module, using the trial tibial locking screw (9066.47.190) (Fig. 25).

The position of the trial module must correspond to the previously determined reference (ex. reference n°6) (Fig. 26).

If the reamer has reached the reference L, the additional L extension must be tightened to the S module (Fig. 26).

Select the appropriate trial stem diameter and tighten it to the module or to the L extension (Fig. 27).

The trial stem diameter corresponds to the one of the last reamer used (see Tab. 3).

Tab. 3

| Stem diameters & length | Stem + S Module 30mm (N, +3, +6) | Stem + S Module 30mm (N, +3, +6) + L extension 25mm (total 55mm) |
|-------------------------|----------------------------------|------------------------------------------------------------------|
| Dia. 14mm length 60mm   | L = 90mm                         | L = 115mm                                                        |
| Dia. 16mm length 60mm   | L = 90mm                         | L = 115mm                                                        |
| Dia. 18mm length 85mm   | L = 115mm                        | L = 140mm                                                        |
| Dia. 20mm length 85mm   | L = 115mm                        | L = 140mm                                                        |
| Dia. 22mm length 110mm  | L = 140mm                        | L = 165mm                                                        |
| Dia. 24mm length 110mm  | L = 140mm                        | L = 165mm                                                        |

If necessary, connect under the surface of the trial tibial plate the trial tibial augment of the corresponding size with the appropriate thickness. The augment is magnetically attached to the trial tibial plate (see page 18 for the preparation).

**NOTE.** To perform a CCK, be aware of using the trial components without the "H" marking. Using the "H" trial tibial plate would lead to misalignment since the hinged system has a 0° tibial slope instead of 6° (CCK).

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## CCK Femur



Figure 28




Figure 29

Figure 29b




Figure 30

### ▼ FEMUR

#### FEMORAL SIZE DETERMINATION

Use the phantoms #1-5 to trial the size of the femur (9066.47.410-450) to approximate femoral size (Fig. 28).

#### FEMORAL STEM SIZE DETERMINATION

Insert the handle on the reamer shaft. Introduce the smallest reamer (14mm) into the femoral canal. Progressively ream the canal until cortical contact is achieved. (Fig. 29).

If the first reference (S) is reached, it means the short trial stem needs to be used together with the correct trial offset module for the femur (+3/-3 L/R) (Fig. 29b) (9066.47.570-595).

If the second reference (L) is reached, it means the long trial stem needs to be used, in combination with the long extension module. Therefore for the final implant assembly, the definitive stem will be assembled with the "L" module (see Tab. 4).

Tab. 4

| Ref | Length | Trial                          | Definitive       |
|-----|--------|--------------------------------|------------------|
| S   | 30mm   | S Femoral Module               | S Femoral Module |
| L   | 55mm   | S Femoral Module + L Extension | L Femoral Module |

In case of sinking over the L reference, proceed with bigger diameter reamers. The last used reamer corresponds to the trial and definitive stem.

**NOTE.** Take note of the diameter and of the depth level (S or L) reached.

In order to avoid any interference between femoral module and bone, ream again the femur using the last reamer until the stop reference line is reached (Fig. 30).

Care should be taken to choose the correct diameter of the stem to avoid pain in the diaphyseal part of the bone.




Figure 31




Figure 32




Figure 33

### DISTAL RESECTION

To re-cut, the distal resection guide (9066.15.050) + the distal re-cut spacer of 9mm (9066.47.480) + the angle guide 6° 9066.47.470 need to be connected over the reamer (Figs. 31-32).

Leaving the reamer in situ, slide onto it the distal resection mask. Place the mask on the existing distal surface. Lock the angled guide by turning the screw.

Secure the jig with pins and proceed with the distal resection by introducing the oscillating saw-blade through the +3mm slot (Fig. 33).

### FEMORAL ROTATION

Check the rotation to position the distal resection guide in line with the epicondylar axis and verify the anatomical landmarks in combination with releases when necessary. The tibial shaft axis can be used as a reference for femoral rotation.

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## CCK Femur




Figure 34

In case of relevant bone loss, it can be difficult to position the distal resection guide. The medial epicondyle can then be used as a reference, by introducing the femoral joint line guide (9066.47.064) into the slot of the distal resection guide (Fig. 34).

If the epicondyle is missing please refer to other anatomical landmarks to guide you in finding the right position.

In order to position the distal resection level at the correct distance from the femoral distal surface, slide the distal resection guide until the joint line guide is touching the medial epicondyle (Fig. 35).

Secure the guide with pins and proceed with the distal resection introducing the saw blade 1.27mm into the +3mm slot (Fig. 33).

Remove the pins and the guide, leaving the reamer into the canal.




Figure 35

### FLEXION GAP & LIGAMENT TENSION

To check the flexion gap after the distal cut, the blue blocks for the trial ligament tension thickness (9066.20.710-750) can be positioned on the tibia. This gives you an idea about the flexion gap, the rotation of the femur and the ligament tension (Fig. 36) when trialing, please remove the reamer & cutting blocks.



Figure 36

### ANTERIOR, POSTERIOR AND CHAMFER RESECTIONS

Take the 4-in-1 femoral guide (9066.47.510-550) of the appropriate size, together with the trial Rt-Lt /Rt+3 Lt -3/ Rt -3 Lt +3 short femoral modules (9066.47.570/575/580/585/590/595) (L or R) and insert both over the reamer defining the correct offset (*Fig. 37*).

With the MULTIGEN-PLUS system the average distance between the center of the medial epicondyle and the distal surface is ~19mm (*Fig. 38*).




Figure 37




Figure 38

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## CCK Femur




Figure 39

In case of relevant bone loss, to aid stability of the 4-in-1 femoral guide, it's possible to insert on the back of the guide two trial distal augments (5mm and 10mm) (9066.47.650/655), guided by the dedicated holes (Figs. 39-40).

Align the 4-in-1 femoral guide to the transepicondylar line or rotate it until the lower part of the femoral guide surface is parallel to the resected posterior condyles.

**NOTE.** Be aware that for size 1, only augment 5mm is available.




Figure 40

### FEMORAL OFFSET DETERMINATION

Introduce a feeler blade or a sickle through the anterior slot of the 4 in 1 guide to evaluate the A/P position of the anterior cut.

This allows to understand if the neutral femoral alignment guide assures the right positioning.

It is possible to move the 4 in 1 guide anteriorly (in order to avoid notching) or posteriorly (getting nearer to the anterior cortex) substituting the neutral femoral alignment guide with an offset femoral alignment guide (+3mm, -3mm) (Figs. 41-42).

Each femoral alignment guide corresponds to one different trial and definitive femoral module (See Tab.5).

Tab. 5

| Femoral alignment guide | Femoral module   |
|-------------------------|------------------|
| R-L                     | R / L Module     |
| R+3 or L-3              | R+3 / L-3 Module |
| R-3 or L+3              | R-3 / L+3 Module |

R = Right; L = Left

Secure the femoral guide and proceed with the anterior, posterior and chamfers resections.




Figure 41




Figure 42




Figure 43

### FEMORAL CCK/H BOX RESECTION

Connect the CCK/H box mask onto the 4-in-1 guide. Secure the screw for CCK/H box built-in the CCK/H box mask (Fig. 43).

To improve stability of the CCK/H box mask in case of bone loss on the anterior part of the femur, insert the screw for CCK/H box (9066.47.071) and tighten it until touching the reamer shaft (Fig. 44).

Insert a narrow saw blade 1.27mm through the slots and proceed with the box preparation (Fig. 45).




Figure 44




Figure 45

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## CCK Femur




Figure 46




Figure 47

**NOTE.** If the distal resection is not possible, use directly the femoral guide.



### **WARNING**

It is important to know exactly the correct position of the femoral guide respect to the medial epicondyle; to define this position insert the joint line guide for femoral box mask (9066.47.065) (R02) into the femoral guide (Fig. 46).

Once the joint line guide is touching the medial epicondyle, secure the 4-in-1 femoral guide with pins and proceed with the preparation of the box introducing the saw blade into the slots (Fig. 47).

Remove the joint line guide for femoral box, the femoral box mask, the 4-in-1 femoral guide and the reamer. Finish the box preparation.



Figure 48

### FEMORAL TRIALS


Connect the CCK/H trial femoral component #1-#5 (9066.20.015/025/035/045/055) of the appropriate size with the femoral module (see Tab. 5) using the femoral locking screw (9066.47.195) (Fig. 48).

*Positioning check:* the mark of the module must be read looking the femoral component from the posterior side (Fig. 49).

If the reamer has reached the reference L, the additional L extension must be tightened to the S module (Fig. 50).

Fasten the trial stem with the trial femoral module (or with the L Extension) (Fig. 51).

Impact the femoral trial components with the provided impactor (Fig. 52).



Left knee (L, L+3, L-3)

Figure 49



Figure 51



Figure 50



Figure 52

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## CCK Femur




Figure 53

### FEMORAL AUGMENTS

Once the femoral trial component is on place, check if it is necessary to perform the resections to host the distal and posterior augments (*Fig. 55*).

Perform the resections (*Figs. 53-54*) and remove the components to assemble the trial distal /posterior femoral augments #1-#5 (9066.71.105/205/210/305/310/405/410/505/510) (9066.72.105/205/210/305/310/405/410/505/510).

Impact the femoral trial component.



Figure 54




Figure 55



Figure 56

### ▼ TRIAL REDUCTION

#### TIBIAL TRIALS WITH LINER

Introduce all the previously assembled trial components into the tibial canal.

**NOTE.** Before inserting the definitive implants, trial reduction is required (Fig. 57).

Insert the CCK cam (9066.35.137) into the tibial trial liner with the appropriate thickness (Fig. 56).

Insert the liner onto the trial tibial plate (Fig. 57).




Figure 57

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## CCK Trial Reduction




Figure 58



Figure 59



Figure 60

### ARTICULAR TEST

Perform the final trial reduction. Check the position of the components and the liner thickness (9066.35.110-610) (Fig. 58).

After the trial reduction remove first the trial tibial liner with the plier for Trial Tibial Liner (9066.35.610) (Fig. 59) and then remove all the components without disassembling them.

Use the multifunction extractor (9066.25.190) to remove the tibial trials assembly (Fig. 60).

**NOTE.** To prevent the tibial plate from shifting the position you can, fix it with 2 long pins for Trial Tibial Plates 9069.10.285 (M9).



Figure 61

### ▼ DEFINITIVE IMPLANT

Open the definitive components (CCK) and assemble them.

#### TIBIAL IMPLANT

The CCK tibial module +3/+6 must be assembled with the definitive cemented fixed tibial plate in the same position determined during the trial.

Overturn the definitive tibial plate attaching the two parts of the tibial module positioning jig (9066.48.190) around the tibial keel (magnetically) (Fig. 61).

Align the CCK tibial module reference mark to the number engraved on the positioning jig which corresponds to the offset position, previously determined during the trial (Fig. 62).

Assemble the CCK tibial module onto the tibial plate by hammering it once, in order to couple the Morse taper. Assemble the appropriate CCK/H stem to the CCK tibial module using the same procedure (Fig. 63).




Figure 62



Figure 63

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### CCK Definitive Implant



Figure 64

**NOTE.** If a tibial augment is utilized, position the augment on the underside of the tibial plate and secure it with two screws from the top using the hexagonal screwdriver 3.5mm (Fig. 64).

Use the screws contained in the package to secure the augment to the tibial plate.

Screw slightly both of them before fully tightening.

Apply a layer of bone cement on the underside of the definitive tibial component.

Carefully insert the definitive CCK tibial component avoiding malrotation. When the implant is fully inserted ensure the fixation by hammering with the tibial impactor (9066.25.110) (Fig. 65).

Remove all the cement residual.




Figure 65



Figure 66

### FEMORAL IMPLANT

Use the femoral module stems coupling base (9066.49.900) to assemble the definitive femoral component with the chosen femoral module.

*Femoral module positioning check:* as previously done with the femoral trials, the mark on the module must be read looking the CCK femoral component from the posterior side (Fig. 66).

Assemble the femoral module onto the CCK femoral component by hammering it once in order to couple the morse taper (Fig. 67).

Then assemble the appropriate CCK/H stem to the femoral module using the same procedure (Fig. 68).



Figure 67



Figure 68

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## CCK Definitive Implant



Figure 69

If distal or posterior augments are used, connect the augments onto the CCK femoral component. Insert the screw through the augments and tighten with the cardan hexagonal screwdriver 3.5 mm (9095.10.223, 9095.10.224) (Fig. 69).

**NOTE.** Use the screw contained in the package to secure the femoral augment to the femoral component.

Place a layer of cement on the underside of the CCK femoral prosthesis

Insert the CCK femoral component onto the distal femur.



### WARNING

Make sure to avoid scratching on the definitive implant.

Be sure that soft tissue is not trapped beneath the implant. Use the Femoral Impactor to fully seat the femoral component (Fig. 70).

Check the medial and lateral sides to make sure the CCK femoral component is fully impacted.

Remove any cement particulate or presence of soft tissue from the definitive tibial plate. Insert onto the tibial plate the appropriate definitive tibial liner.

First, slide the definitive CCK liner posteriorly onto the tibial plate in order to insert the polyethylene back lip beneath the posterior tooth of the tibial plate.

Then connect the liner anteriorly with continuous pressure to avoid a rebound of the cemented implant.

As an alternative wait until the cement has hardened to impact the liner with the impactor (9066.30.160) (Fig. 71).



Figure 70



Figure 71

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### CCK Definitive Implant



Figure 72

**NOTE.** Position the knee in light flexion (10°) to let harden the cement out (Fig. 73)

A locking screw is required for the CCK tibial liner. Insert the screw through the tibial liner post and then strongly tighten it using the hexagonal screwdriver 3.5mm (Fig. 72).

Proceed with the CCK definitive implant (Fig. 73) move the yellow part to the previous part: and position the knee in light flexion (10°) to let harden the cement out.

**NOTE.** The locking screw is packaged with the CCK tibial liner. The locking screw has an anti-unscrewing system which keeps the screw itself completely inside the tibial liner post.



Figure 73





Figure 1

### ▼ TIBIA

#### TIBIAL STEM SIZE DETERMINATION

Introduce the first smallest reamer (14mm) into the tibial canal (Figs. 1-2).

Progressively ream the canal until cortical contact is achieved.

If the stability is reached at the reference number S, it means that a short stem is needed (Fig. 3).

The trial short stem (9066.47.230 - 9066.47.255) needs to be assembled in combination with the correct tibial offset trial module (0/+3/+6).

For the definitive implant you should take the correct stem diameter + SHORT tibial module straight/+3mm/+6mm.



Figure 2




Figure 3

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### H Tibia




Figure 4

If the second line (L) is reached, it means a long stem is needed (*Fig. 4*). Therefore use the trial tibial stem connected with the chosen tibial trial offset module (straight/+3mm/+6mm and the extension module (long) (9066.47.615). The final implant will consist of the definitive stem and the “L” module (*Tab. 1*).

Tab. 1

| Ref      | LENGTH      | TRIAL                            | DEFINITIVE      |
|----------|-------------|----------------------------------|-----------------|
| <b>S</b> | <b>30mm</b> | S Tibial Module                  | S Tibial Module |
| <b>L</b> | <b>55mm</b> | S Tibial Module +<br>L Extension | L Tibial Module |

Whether the reamer goes below the L reference use reamers with the highest reference point as diameters.

If there is a difference between the medial and the lateral plateau use the highest as reference.

The size of the latest reamer used corresponds to the diameter of both trial and definitive stem.

When between 2 sizes please chose the smaller one to avoid pain in the diaphyseal part of the bone.




Figure 5

### TIBIAL RESECTION (HINGE 0° OF SLOPE)

Remove the handle from the reamer and insert the prismatic guide onto the reamer shaft.

Insert the IM tibial cutting guide H (9066.48.050) on the prismatic guide.

**NOTE.** The H tibial resection jig has a 0° tibial slope.

Introduce the sickle (9066.12.010) or a free blade through the cutting guide slot and place it on the highest tibial plateau (if there is any difference between the medial and the lateral plateau) (Fig. 5).

Lock the screw of the prismatic guide on the reamer arm. Rotate the millimetric screw until the desired resection level is reached and secure the jig using pins.

It is also possible to use the +0mm tibial depth resection (H2) (9066.25.160) gauge to measure the gap.

Proceed with the resection (Fig. 6).




Figure 6

Whether a tibial augment is needed it's possible to directly perform the +7mm or the +12mm resection to host it. Rotate the millimetric screw until the desired resection level +7mm or +12mm, referring to the last performed resection, is reached (Fig. 7).

Remove the pins and the resection jig.

If needed, it is possible to perform the tibial augments cuts after the tibial resection mask has been removed, using the augment tibial resection guide that has to be inserted onto the reamer shaft as well (Fig. 12) on page 18.




Figure 7



optional  
9066.48.055



2 different trial tibial augments can be used:  
7mm (9066.47.315/325/335/345/355)  
12mm (9066.47.320/330/340/350/360)

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### H Tibia




Figure 8

#### TIBIAL PLATE SIZE AND TIBIAL STEM OFFSET DETERMINATION

Leaving the reamer in situ, place the trial tibial plate (9066.48.110-9066.48.150) on the tibial surface.

Slide the straight tibial dialer (9066.47.170) along the shaft and insert it into the circular slot of the trial tibial plate.

At the straight tibial dialer corresponds the trial short straight tibial module (without offset) (9066.47.200).

Select the tibial size to achieve maximal tibial coverage (Figs. 8-9).

**NOTE.** In a Hinge surgery be aware of using the tibial trial plates with marking H.




Figure 9




Figure 10

If a tibial augment is required, a trial augment (9066.47.315-360) of the same size of the tibial plate can be fixed magnetically to the trial plate.

To obtain a better coverage of the resected tibial surface, it could be necessary to offset the trial tibia plate replacing the neutral dialer with the +3mm or +6mm eccentrical dialer and rotating it until the best positioning is reached (Fig. 10).

A +3mm and +6mm dialer corresponds to a +3mm and +6mm tibial module (see Tab. 2).

Tab. 2

| Dialer | Trial/Definitive                      |
|--------|---------------------------------------|
| 0      | trial short straight tibial module    |
| +3     | trial short offset +3mm tibial module |
| +6     | trial short offset +6mm tibial module |

**NOTE.** Once the right positioning is reached, take note of the number marked on the dialer aligned with the reference on the trial tibial plate (Fig. 11).




Figure 11

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## H Tibia



Figure 12




Figure 13




Figure 14

### TIBIAL KEEL PREPARATION

Insert two K-wires (dia. 2mm) through the small diameter holes of the H trial tibial plate (Fig. 12).

**NOTE.** The K-wires are not included in the set.

Remove the dialer, the H trial tibial plate and the reamer.

Reposition the H trial tibial plate using the two K-wires as guides (Fig. 13) and secure the H tibial plate by using the pins.

If a H trial tibial augment is used, secure the plate using the long pins.

Remove the two K-wires.


Attach the guide for the “H” Tibial Winged Broach guide (9066.47.095) onto the H trial tibial plate (Fig. 14).

**NOTE.** Make sure to use the winged broach for Hinge, the one with the H mark.

To check the tibial alignment, the alignment rod (9066.15.090) (Fig. 15) can be used through the slot of the handle.

Figure 15





Insert the tibial winged broach (9066.30.100) into the "H" Tibial winged broach guide (9066.47.095) until the lever tooth engages (Fig. 16).

Use the multifunction extractor (9066.25.190) to remove the "H" tibial winged broach guide and the guide (Fig. 17).

In order to avoid any interference between tibial module and bone, reattach the guide for the broach onto the trial tibial plate and ream the tibia using the tibial reamer (9066.47.225) for until it's in even contact with the guide for the tibial broach (Fig. 18).

This procedure is necessary when a tibial module with +6mm offset needs to be used and with +3mm offset tibial module when the last reamer is smaller than 20mm.




Figure 16



Tibial reamer with stop  
9066.47.225



Figure 18

Figure 17

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## H Tibia



Figure 19




Figure 20



Figure 21

### TIBIAL TRIALS

Connect the H trial tibial plate with the H trial tibial keel and the chosen tibial module, using the tibial locking screw (Fig. 19).

The position of the trial module must correspond to the previously determined reference (ex. reference n°6) (Fig. 20).

If the reamer has reached the reference L, the additional L extension must be tightened to the S module (Fig. 20). Select the appropriate trial stem diameter and tighten it to the module or to the L extension (Fig. 21).

The trial stem corresponds to the last reamer used (see Tab. 3).

Tab. 3

| Stem diameters & length | Stem + S Module 30mm (N, +3, +6) | Stem + S Module 30mm (N, +3, +6) + L extension 25mm (total 55mm) |
|-------------------------|----------------------------------|------------------------------------------------------------------|
| Dia. 14mm length 60mm   | L = 90mm                         | L = 115mm                                                        |
| Dia. 16mm length 60mm   | L = 90mm                         | L = 115mm                                                        |
| Dia. 18mm length 85mm   | L = 115mm                        | L = 140mm                                                        |
| Dia. 20mm length 85mm   | L = 115mm                        | L = 140mm                                                        |
| Dia. 22mm length 110mm  | L = 140mm                        | L = 165mm                                                        |
| Dia. 24mm length 110mm  | L = 140mm                        | L = 165mm                                                        |

If necessary, connect under the surface of the H trial tibial plate the augments of the corresponding size with the appropriate thickness. The augment is magnetically attached to the H trial tibial plate.

**NOTE.** To perform a H, be aware of using the trial components with the “H” marking.

Using the CCK trial tibial plate would lead to misalignment since the CCK system has a 6° tibial slope instead of a 0° for Hinge.



Figure 22

### ▼ FEMUR

#### FEMORAL SIZE DETERMINATION

Use the lateral femoral templates (9066.47.410-450) to approximate femoral size (Fig. 22).

#### FEMORAL STEM SIZE DETERMINATION

Insert the handle on the reamer shaft. Introduce the smallest reamer (14mm) into the femoral canal. Progressively ream the canal until cortical contact is achieved. (Fig. 23).




Figure 23

Figure 23b

If the first reference (S) is reached, it means the short trial stem needs to be used together with the correct trial offset module for the femur (+3/-3 L/R) (Fig. 23b).

If the second reference (L) is reached, it means the long trial stem needs to be used, in combination with the long extension module. Therefore for the final implant assembly, the definitive stem will be assembled with the "L" module (see Tab. 4).

Tab. 4

| Ref | Length | Trial                          | Definitive       |
|-----|--------|--------------------------------|------------------|
| S   | 30mm   | S Femoral Module               | S Femoral Module |
| L   | 55mm   | S Femoral Module + L Extension | L Femoral Module |

In case of sinking over the L reference, proceed with bigger diameter reamers.

The last used reamer corresponds to the trial and definitive stem.

**NOTE.** Take note of the diameter and of the depth level (S or L) reached.

In order to avoid any interference between femoral module and bone, ream again the femur using the last reamer until the stop reference line is reached (Fig. 24).




Figure 24

When between 2 sizes please chose the smallest one to avoid pain in the diaphyseal part of the bone.

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## H Femur




Figure 25




Figure 26




Figure 27

### DISTAL RESECTION

To re-cut the distal surface, insert on the distal resection guide the re-cut distal spacer (9066.47.480) (*Fig. 25*).

The distal resection guide has a 6° of valgus angle build in (*Fig. 26*) and is used for a left and a right knee. M/L positioning is done by the symmetrical sizing.

Leaving the reamer in situ, slide onto it the distal resection guide (9066.15.050) with the angled guide 6° (9066.47.470) in the correct position for a left or a right knee and position on the existing distal surface. Lock the angled guide by turning the screw.

Secure the distal resection guide with pins and proceed with the distal resection by introducing the oscillating saw-blade through the +3mm slot (*Fig. 27*).




Figure 28

In case of relevant bone loss, it can be difficult to position the distal resection mask. The medial epicondyle can then be used as a reference, by introducing the joint line guide for femoral distal mask (9066.47.064) into the slot of the distal resection guide (9066.15.050) (Fig. 28).

In order to position the distal resection level at the correct distance from the femoral distal surface, slide the distal resection guide until the joint line guide is touching the medial epicondyle (Fig. 29).

Secure the guide with pins and proceed with the distal resection introducing the saw blade into the +3mm slot (Fig. 27).

Remove the pins and the jig, leaving the reamer into the canal.




Figure 29

### FLEXION GAP & LIGAMENT TENSION

To check the flexion gap after the distal cut, the blue blocks for the trial ligament tension thickness (9066.20.710-750) can be positioned on the tibia. This gives you an idea about the flexion gap, the rotation of the femur and the ligament tension (Fig. 30).

Please remove the reamer and cutting blocks when trialing.



Figure 30

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## H Femur




Figure 31

### ANTERIOR, POSTERIOR AND CHAMFER RESECTIONS

Take the femoral guide (9066.47.510-550) of the appropriate size, together with the neutral femoral alignment guide (L or R) (9066.47.570-9066.47.585) and insert it over the reamer (*Fig. 31*).

With the MULTIGEN-PLUS system the average distance between the center of the medial epicondyle and the distal surface is ~19 mm (*Fig. 32*).




Figure 32




Figure 33

In case of relevant bone loss, to aid stability of the femoral guide, it's possible to insert on the back of the mask two trial distal augments (5mm and 10mm) (9066.47.650/655), guided by the dedicated holes (Fig. 33-34).

Align the femoral guide to the transepicondylar line or rotate it until the lower mask surface is parallel to the resected posterior condyles.

**NOTE.** Be aware that for size 1, only augment 5mm is available.




Figure 34

### FEMORAL OFFSET DETERMINATION

Introduce a feeler blade or a sickle through the anterior slot of the femoral guide to evaluate the A/P position to the anterior cut.

This allows to understand if the neutral femoral alignment guide assures the right positioning.

It is possible to move the femoral guide anteriorly (in order to avoid notching) or posteriorly (getting nearer to the anterior cortex) substituting the neutral femoral alignment guide with an offset femoral alignment guide (+3mm, -3mm) (9066.47.575/580/590/595) (Fig. 35-36).

Each femoral alignment guide corresponds to one different trial and definitive femoral module (See Tab.5).

Tab. 5

| Femoral alignment guide | Femoral module   |
|-------------------------|------------------|
| R-L                     | R / L Module     |
| R+3 or L-3              | R+3 / L-3 Module |
| R-3 or L+3              | R-3 / L+3 Module |

R = Right; L = Left

Secure the femoral guide and proceed with the anterior, posterior and chamfer resections.




Figure 35





Figure 36

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## H Femur





### FEMORAL CCK/H BOX RESECTION

Connect the CCK/H box mask (9066.47.070) onto the femoral guide.

Secure the CCK/H box with the screw built-in the femoral box mask (*Fig. 37*).

To improve stability of the CCK/H box mask in case of bone loss on the anterior part of the femur, insert the screw for CCK/H box (9066.47.071) and tighten it until touching the reamer shaft (*Fig. 38*).

Insert a narrow saw blade through the slots and proceed with the box preparation (*Fig. 39*).






Figure 40

**NOTE.** If the distal resection is not possible, use directly the femoral guide.



### WARNING

It is important to know exactly the correct position of the femoral guide in respect to the medial epicondyle; to define this position, insert the joint line guide for femoral box mask (9066.47.065) into the femoral guide (Fig. 40).

Once the joint line guide is touching the medial epicondyle, secure the femoral guide with pins and proceed with the preparation of the box introducing the saw blade into the slots (Fig. 41).

Remove the joint line guide for femoral box, the femoral box mask, the femoral guide and the reamer. Finish the box preparation.





Figure 41

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## H Femur



Figure 42



Left knee (L, L+3, L-3)

Figure 43



Figure 44

### FEMORAL TRIALS

Connect the trial femoral component (9066.20.015-055) of the appropriate size with the trial femoral module (see Tab. 5) using the femoral locking screw ((9066.47.195) (Fig. 42).

*Positioning check:* the mark of the module must be read looking the femoral component from the posterior side (Fig. 43).

If the reamer has reached the reference L, the additional L extension ((9066.47.615) must be tightened to the S module (Fig. 44).

Fasten the trial stem with the trial femoral module (or with the L Extension) (Fig. 45).

Impact the femoral trial components with the provided impactor (9069.10.220).



Figure 45



Figure 46




Figure 47

### FEMORAL AUGMENTS

Once the femoral trial component is on place, check if it is necessary to perform the resections to host the distal and posterior augments (*Fig. 47*).

Perform the resections (*Figs. 48-49*) and remove the components to assemble the augments.




Figure 48




Figure 49

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## H Trial Reduction




Figure 50




Figure 51

### ▼ H TRIAL REDUCTION

#### TIBIAL TRIALS WITH LINER

Introduce all the previously assembled trial components into the tibial canal, and gently impact them with the tibial impactor (9066.25.110).

Insert the CCK/H tibial trial liner with the appropriate thickness and balance the gap (*Fig. 50-51*).

**NOTE.** To prevent translational forces, tilting and subluxation with the trials you can use the CCK peg to trial with the trial femur as well. In case of trialing with the definitive Hinge implant, you can't use the peg because of the Hinge mechanism.



Figure 52

### ARTICULAR TEST

Perform the final trial reduction. Check the position of the components and the liner thickness (*Fig. 52*).

After the trial reduction, remove all the components without disassembling them.

Use the multifunction extractor to remove the tibial trials assembly (*Fig. 53*).

**NOTE.** For Hinged you can use the definitive femoral component in combination with the trial liner without cam.

For Hinged you can also use the definitive femoral component in combination with the trial liner without cam.



Figure 53

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### H Definitive implant



Figure 54

#### H TIBIAL IMPLANT

The H tibial module with offset (straight/+3/+6) must be assembled with the definitive H tibial plate in the same position determined during the trial.

Use the definitive H tibial plate and attach the two parts of the tibial module positioning jig (9066.48.190) around the tibial keel, (the 2 parts are magnetic) (Fig. 54).

Align the H tibial module with reference mark to the number engraved on the positioning jig, which corresponds to the offset position previously determined during the tibial preparation (Fig. 55).




Figure 55

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## H Definitive implant



Figure 56

Assemble the H tibial module onto the H tibial plate by striking it until the connection is stable with a hammer, in order to couple the Morse taper.

Using the same procedure, assemble the appropriate CCK/H stem to the H tibial module (Fig. 56).

If a H tibial augment is utilized, firmly fasten the screws to the tibial plate (Fig. 57) using the following procedure:

1. the screws must be completely tightened to the H tibial augment;
2. the protruding part of the screws must be inserted into their seat on the H tibial plate, and then the H tibial augment must be pressed on the H tibial plate.

**NOTE.** Use the two screws contained in the package to secure the H tibial augment to the tibial plate.



Figure 57

Apply a layer of bone cement to the underside of the definitive H tibial component.

Carefully insert the definitive H tibial component avoiding malrotation. When fully inserted, several hammer blows might be delivered to the end of the tibial impactor. Remove all residual cement (Fig. 58).




Figure 58

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## H Definitive implant



Figure 59



Figure 60



Figure 61

### H FEMORAL IMPLANT

Use the femoral module stems coupling base (9066.49.900) to assemble the definitive appropriate sized H femoral component with the appropriate femoral module.

*Femoral module positioning check:* as previously done with the femoral trials, the mark of the module must be read looking the H femoral component from the posterior side.

Assemble the femoral module onto the H femoral component by hammering once in order to couple the morse taper (*Fig. 59*).

Using the same procedure, assemble the appropriate CCK/H stem to the femoral module (*Fig. 60*).

If distal or posterior augments are used, attach the definitive augment onto the H femoral component. Insert the screw through the augment and tighten with the Hexagonal Screwdriver 3.5 mm (cardan joint) (*Fig. 61*).

**NOTE.** Use the screw contained in the package to secure the femoral augment to the H femoral component.

Place a layer of cement on the innerside of the H femoral prosthesis. Insert the H femoral component onto the distal femur. Make sure to avoid scratching the implant component surfaces.

Be sure that soft tissue is not trapped beneath the implant. Use the Femoral Impactor to fully seat the H femoral component (*Fig. 62*). Check the medial and lateral sides to make sure the H femoral component is fully impacted.




Figure 62




Figure 63

### FEMORAL BOX CHECK

Before opening the definitive implant, insert the femoral box gauge (9066.48.220) between the femoral condyles, and check if the central distal resection is deep enough. Then, place the definitive implant following the instructions below (Fig. 63).

### ▼ H HINGE FINAL ASSEMBLY



#### WARNING

It is very important to follow the procedure step by step.



#### WARNING

The use of the code 9066.51.000 (CCK-H System Additional Set - Multigen Plus) in combination with the code 9095.11.751 (Torque wrench) IS MANDATORY.

Open the definitive H tibial liner.

Take the polyethylene component only, and place it onto the metal peg on the H tibial plate (Fig. 64).

Insert the H Hinge through the H tibial liner into the hole on the H tibial plate (Fig. 65).

While distracting the joint, flex the femoral Hinge Post to align it perfectly with the Morse taper of the H Hinge (Fig. 66).



Figure 64



Figure 65



Figure 66

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### H Definitive implant




Figure 67

#### PREPARATION PHASE

Insert the threaded shaft of the tibial pin positioner (9066.48.205) through the hole of the hinge post and screw it gently to the H-Hinge until it stops automatically (Fig. 67).

**NOTE.** "Finger tight" do not over tighten.

Insert the femoral counter Tourque (9066.48.200) on the threaded shaft and flush it against the femoral condyles, holding it firmly (Fig. 68).




Figure 68

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## H Definitive implant




Figure 69

Slide the ferrule of the tibial pin positioner second part of the tibial pin positioner (9066.48.205) the threaded shaft and screw it clockwise until it is locked (*Fig. 69*).

This procedure will move the Morse cone upwards into the female cone of the hinge coupling the two components (*Fig. 70*).

Insert the tibial pin positioner adapter (9066.48.210) (*Fig. 71*) along the threaded shaft of the tibial pin positioner and couple its hex extremity with the hexagonal slot of the upper part of the ferrule.



### WARNING

Do not push downwards during this step.




Figure 70




Figure 71

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### H Definitive implant



Figure 72

Place the torque wrench (9095.11.751) on the tibial pin positioner adapter (*Fig. 72*).

**NOTE.** The torque wrench 9095.11.751 is not included for quality reasons in the set 9066.51.000.

Turn the torque wrench clockwise until a “click” is heard (*Fig. 73*).

Loosen the torque wrench a half-turn counterclockwise to remove it from the tibial pin positioner adapter (*Fig. 74*). Remove the torque wrench and the tibial pin positioner adapter.



Figure 73



Figure 74




Figure 75

Remove the tibial pin positioner ferrule by turning it counterclockwise (*Fig. 75*).

In order to remove the threaded shaft, use the hexagonal slot of the pin positioner adapter to engage its extremity (*Fig. 76*).

Then unscrew the threaded shaft using the adapter as a lever. Eventually remove all the instruments (*Fig. 77*).




Figure 76




Figure 77

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### H Definitive implant



Figure 78



Figure 79



Figure 80

#### FINAL PHASE

Use the dynamometric screwdriver 3.5mm 5Nm (9095.10.226) to insert the H locking screw through the hinge post and then tighten it until the 5 mark reaches the reference line (Fig. 78).

Remove the dynamometric screwdriver (Fig. 80).

**NOTE.** 5Nm is the required torque in order to fully tighten the screw, ensuring the system to be safe, well assembled and locked. This value is marked on the dynamometric screwdriver: it is recommended not to exceed this value to avoid any damage on the screw (Fig. 79).

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Extraction of Hinge components



Figure 81

### EXTRACTION

To extract the implant flex the knee and start to remove the top screw by using the dynamometric screwdriver (*Fig. 81*).


Insert the threaded shaft of the tibial pin positioner through the hinge, into the hinge post and thread it ("finger tight" do not over tighten) (*Fig. 82*).



Figure 82

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### Extraction of Hinge components



Slightly hammer on it to disconnect the taper (*Fig. 83*).

Extract the threaded shaft and uncouple the parts to be removed (*Fig. 80*).

Figure 83




Figure 84




Figure 1

### ▼ PATELLAR PROSTHESIS

#### MEASUREMENT OF THE PATELLA

Using the patellar gauge, measure the thickness of the patella before the resection (*Fig. 1*).

#### PATELLAR RESECTION

If the patella is 20mm thick at least, let the stylus marked 10mm STD rest on the patella. In this case the amount of bone resection will be 10mm, the same thickness as the definitive patellar prosthesis (*Fig. 2*).

If the measured thickness is lower than 20mm, in order to maintain at least 10mm of bone, the amount of bone resection must be lower than 10mm.

Let the 8mm or 6mm stylus rest on the patella. In this case the amount of bone resection and the thickness of the definitive patellar prosthesis will be 8 or 6mm.

For the resection, insert the saw blade through the slot from the lateral side of the plier (*Fig. 3*).




Figure 2



Figure 3

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### Patellar prosthesis




Figure 4




Figure 5



Figure 6

#### FINAL PATELLA PREPARATION

Place and align the appropriate sized patellar drilling guide to the center of resected patellar surface and secure it (*Fig. 4*).

Position the patellar drilling guide on the patellar mask and drill the peg holes using the patellar drill (*Figs. 5-6*).


Place the trial patellar component onto the resected patellar surface and check the patellar tracking.

**NOTE.** A domed symmetric patella should always be placed on the medial border of the bony cortex to avoid tilting. The full coverage is not possible and needed, while a proper alignment is mandatory.


# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Instrument Set

### ▼ 9066.47.000 MULTIGEN-PLUS CCK-H Common Set n. 8



Upper tray




Lower tray

| Ref. | CODE        | DESCRIPTION                              | Qty. |
|------|-------------|------------------------------------------|------|
| A47  | 9066.47.005 | Handle for Reamers                       | 2    |
| B47  | 9066.47.014 | Reamer Dia. 14mm                         | 1    |
| B47  | 9066.47.016 | Reamer Dia. 16mm                         | 1    |
| B47  | 9066.47.018 | Reamer Dia. 18mm                         | 1    |
| B47  | 9066.47.020 | Reamer Dia. 20mm                         | 1    |
| B47  | 9066.47.022 | Reamer Dia. 22mm                         | 1    |
| B47  | 9066.47.024 | Reamer Dia. 24mm                         | 1    |
| C47  | 9066.47.080 | Unlock Pivot                             | 2    |
| D47  | 9066.47.085 | "L" Hexagonal Wrench 8mm - 3mm           | 1    |
| E47  | 9066.47.090 | M8-M6 Adaptor                            | 1    |
| F47  | 9066.47.095 | "H" Tibial Winged Broach                 | 1    |
| G47  | 9066.47.230 | Trial Stem Dia. 14mm                     | 2    |
| G47  | 9066.47.235 | Trial Stem Dia. 16mm                     | 2    |
| G47  | 9066.47.240 | Trial Stem Dia. 18mm                     | 2    |
| G47  | 9066.47.245 | Trial Stem Dia. 20mm                     | 2    |
| G47  | 9066.47.250 | Trial Stem Dia. 22mm                     | 2    |
| G47  | 9066.47.255 | Trial Stem Dia. 24mm                     | 2    |
| H47  | 9066.47.615 | Extension                                | 2    |
| I47  | 9095.10.223 | Hexagonal Screwdriver 3.5mm              | 1    |
| J47  | 9095.10.224 | Hexagonal Screwdriver 3.5mm Cardan Joint | 1    |
|      | 9066.47.950 | Sterilizable Box                         | 1    |

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Instrument Set


### ▼ 9066.48.000 MULTIGEN-PLUS CCK-H Tibia Set n. 9



| Ref. | CODE        | DESCRIPTION                           | Qty. |
|------|-------------|---------------------------------------|------|
| A9   | 9066.35.137 | CCK-H Cam #1-5                        | 2    |
| B9   | 9066.47.050 | IM Tibial Cutting Guide               | 1    |
| C9   | 9066.47.110 | CCK Trial Tibial Plate #1             | 1    |
| C9   | 9066.47.120 | CCK Trial Tibial Plate #2             | 1    |
| C9   | 9066.47.130 | CCK Trial Tibial Plate #3             | 1    |
| C9   | 9066.47.140 | CCK Trial Tibial Plate #4             | 1    |
| C9   | 9066.47.150 | CCK Trial Tibial Plate #5             | 1    |
| D9   | 9066.47.170 | Straight Tibial Dialer                | 1    |
| D9   | 9066.47.173 | 3mm Eccentric Tibial Dialer           | 1    |
| D9   | 9066.47.176 | 6mm Eccentric Tibial Dialer           | 1    |
| E9   | 9066.47.180 | CCK Tibial Keel                       | 1    |
| F9   | 9066.47.190 | Trial Tibial Locking Screw            | 1    |
| G9   | 9066.47.200 | Trial Short Straight Tibial Module    | 1    |
| G9   | 9066.47.203 | Trial Short Offset +3mm Tibial Module | 1    |
| G9   | 9066.47.206 | Trial Short Offset +6mm Tibial Module | 1    |
| H9   | 9066.47.315 | Trial Tibial Augment #1 h=7mm         | 2    |
| H9   | 9066.47.320 | Trial Tibial Augment #1 h=12mm        | 2    |
| H9   | 9066.47.325 | Trial Tibial Augment #2 h=7mm         | 2    |
| H9   | 9066.47.330 | Trial Tibial Augment #2 h=12mm        | 2    |
| H9   | 9066.47.335 | Trial Tibial Augment #3 h=7mm         | 2    |
| H9   | 9066.47.340 | Trial Tibial Augment #3 h=12mm        | 2    |
| H9   | 9066.47.345 | Trial Tibial Augment #4 h=7mm         | 2    |
| H9   | 9066.47.350 | Trial Tibial Augment #4 h=12mm        | 2    |
| H9   | 9066.47.355 | Trial Tibial Augment #5 h=7mm         | 2    |
| H9   | 9066.47.360 | Trial Tibial Augment #5 h=12mm        | 2    |
| I9   | 9066.48.050 | IM Tibial Cutting Guide H             | 1    |
| J9   | 9066.48.110 | H Trial Tibial Plate #1               | 1    |
| J9   | 9066.48.120 | H Trial Tibial Plate #2               | 1    |
| J9   | 9066.48.130 | H Trial Tibial Plate #3               | 1    |
| J9   | 9066.48.140 | H Trial Tibial Plate #4               | 1    |
| J9   | 9066.48.150 | H Trial Tibial Plate #5               | 1    |
| K9   | 9066.48.180 | H Tibial Keel                         | 1    |
| L9   | 9066.48.190 | Tibial Module Positioning Jig         | 1    |
| M9   | 9069.10.285 | Long Pins for Trial Tibial Plates     | 6    |
|      | 9066.48.950 | Sterilizable Box                      | 1    |

### ▼ 99066.49.000 MULTIGEN-PLUS CCK-H Femur Set n. 10


| Ref. | CODE        | DESCRIPTION                               | Qty. |
|------|-------------|-------------------------------------------|------|
| D10  | 9066.47.071 | Screw for CCK-H box                       | 1    |
| E10  | 9066.47.195 | Trial Femoral Locking Screw               | 2    |
| F10  | 9066.47.410 | Phantom #1                                | 1    |
| F10  | 9066.47.420 | Phantom #2                                | 1    |
| F10  | 9066.47.430 | Phantom #3                                | 1    |
| F10  | 9066.47.440 | Phantom #4                                | 1    |
| F10  | 9066.47.450 | Phantom #5                                | 1    |
| J10  | 9066.47.570 | Right Femoral Alignment Guide             | 1    |
| J10  | 9066.47.575 | Right Offset +3mm Fem. Alignment Guide    | 1    |
| J10  | 9066.47.580 | Right Offset -3mm Fem. Alignment Guide    | 1    |
| J10  | 9066.47.585 | Left Femoral Alignment Guide              | 1    |
| J10  | 9066.47.590 | Left Offset +3mm Fem. Alignment Guide     | 1    |
| J10  | 9066.47.595 | Left Offset -3mm Fem. Alignment Guide     | 1    |
| N10  | 9066.71.105 | Trial Distal Femoral Augment #1 h=5mm     | 2    |
| N10  | 9066.71.205 | Trial Distal Femoral Augment #2 h=5mm     | 2    |
| N10  | 9066.71.210 | Trial Distal Femoral Augment #2 h=10mm    | 2    |
| N10  | 9066.71.305 | Trial Distal Femoral Augment #3 h=5mm     | 2    |
| N10  | 9066.71.310 | Trial Distal Femoral Augment #3 h=10mm    | 2    |
| N10  | 9066.71.405 | Trial Distal Femoral Augment #4 h=5mm     | 2    |
| N10  | 9066.71.410 | Trial Distal Femoral Augment #4 h=10mm    | 2    |
| N10  | 9066.71.505 | Trial Distal Femoral Augment #5 h=5mm     | 2    |
| N10  | 9066.71.510 | Trial Distal Femoral Augment #5 h=10mm    | 2    |
| O10  | 9066.72.105 | Trial Posterior Femoral Augment #1 h=5mm  | 2    |
| O10  | 9066.72.205 | Trial Posterior Femoral Augment #2 h=5mm  | 2    |
| O10  | 9066.72.210 | Trial Posterior Femoral Augment #2 h=10mm | 2    |
| O10  | 9066.72.305 | Trial Posterior Femoral Augment #3 h=5mm  | 2    |
| O10  | 9066.72.310 | Trial Posterior Femoral Augment #3 h=10mm | 2    |
| O10  | 9066.72.405 | Trial Posterior Femoral Augment #4 h=5mm  | 2    |
| O10  | 9066.72.410 | Trial Posterior Femoral Augment #4 h=10mm | 2    |
| O10  | 9066.72.505 | Trial Posterior Femoral Augment #5 h=5mm  | 2    |
| O10  | 9066.72.510 | Trial Posterior Femoral Augment #5 h=10mm | 2    |
| P10  | 6622.15.020 | Femoral Augment Locking Screw             | 5    |



Upper tray

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Instrument set



Lower tray


| Ref. | CODE        | DESCRIPTION                          | Qty. |
|------|-------------|--------------------------------------|------|
| A10  | 9066.20.015 | CCK/H Trial Femoral Component - #1   | 1    |
| A10  | 9066.20.025 | CCK/H Trial Femoral Component - #2   | 1    |
| A10  | 9066.20.035 | CCK/H Trial Femoral Component - #3   | 1    |
| A10  | 9066.20.045 | CCK/H Trial Femoral Component - #4   | 1    |
| A10  | 9066.20.055 | CCK/H Trial Femoral Component - #5   | 1    |
| B10  | 9066.47.060 | Thickness Gauge                      | 1    |
| C10  | 9066.47.070 | CCK/H Box Mask                       | 1    |
| G10  | 9066.47.470 | Angled Guide 6°                      | 1    |
| H10  | 9066.47.480 | Distal Re-Cut Spacer                 | 1    |
| K10  | 9066.47.600 | Trial Rt - Lt Short Femoral Module   | 1    |
| K10  | 9066.47.605 | Trial Rt+3 Lt-3 Short Femoral Module | 1    |
| K10  | 9066.47.610 | Trial Rt-3 Lt+3 Short Femoral Module | 1    |
| I10  | 9066.47.510 | Femoral Guide #1                     | 1    |
| I10  | 9066.47.520 | Femoral Guide #2                     | 1    |
| I10  | 9066.47.530 | Femoral Guide #3                     | 1    |
| I10  | 9066.47.540 | Femoral Guide #4                     | 1    |
| I10  | 9066.47.550 | Femoral Guide #5                     | 1    |
| L10  | 9066.49.900 | Femoral Module Stems Coupling Base   | 1    |
| M10  | 9066.49.910 | Tibial Module Stems Coupling Base    | 1    |
|      | 9066.49.950 | Sterilizable Box                     | 1    |

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Instrument set

### ▼ 9066.55.000 CCK-H Base Set n. 1


| Ref. | CODE        | DESCRIPTION                        | Qty. |
|------|-------------|------------------------------------|------|
| A1   | 9066.15.092 | Dia.3x60mm Fiches                  | 4    |
| A1   | 9066.15.095 | Dia.3x80mm Fiches                  | 4    |
| B1   | 9066.15.235 | Fiches beater                      | 1    |
| D1   | 9066.22.010 | Pin Driver                         | 1    |
| E1   | 9066.22.060 | Screwed Pin with Head Dia.3x60mm   | 4    |
| E1   | 9066.22.080 | Screwed Pin with Head Dia.3x80mm   | 4    |
| F1   | 9066.22.160 | pre-drill Dia. 3mm                 | 1    |
| G1   | 9066.22.170 | Zimmer Rapid Connector             | 1    |
| H1   | 9066.22.180 | Tibial Nail Extractor Plier        | 1    |
| I1   | 9066.22.260 | Pin with Head dia.3x60mm           | 4    |
| I1   | 9066.22.280 | Pin with Head dia.3x80mm           | 4    |
| J1   | 9066.24.060 | Screwed Pin Without Head Dia. 3x60 | 4    |
| J1   | 9066.24.080 | Screwed Pin Without Head Dia. 3x80 | 4    |
| K1   | 9066.30.160 | Liners Impactor                    | 1    |
| M1   | 9066.35.600 | Handle for Trial Tibial Plate      | 1    |



Upper tray

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Instrument set




Lower tray

| Ref. | CODE        | DESCRIPTION                             | Qty. |
|------|-------------|-----------------------------------------|------|
| C1   | 9066.20.710 | h 10mm Trial Ligament Tension Thickness | 1    |
| C1   | 9066.20.720 | h 12mm Trial Ligament Tension Thickness | 1    |
| C1   | 9066.20.730 | h 14mm Trial Ligament Tension Thickness | 1    |
| C1   | 9066.20.740 | h 17mm Trial Ligament Tension Thickness | 1    |
| C1   | 9066.20.750 | h 20mm Trial Ligament Tension Thickness | 1    |
| L1   | 9066.35.110 | CR- PS Trial Liner – #1 / 10mm          | 1    |
| L1   | 9066.35.112 | CR- PS Trial Liner – #1 / 12mm          | 1    |
| L1   | 9066.35.114 | CR- PS Trial Liner – #1 / 14mm          | 1    |
| L1   | 9066.35.117 | CR- PS Trial Liner – #1 / 17mm          | 1    |
| L1   | 9066.35.120 | CR- PS Trial Liner – #1 / 20mm          | 1    |
| L1   | 9066.35.210 | CR- PS Trial Liner – #2 / 10mm          | 1    |
| L1   | 9066.35.212 | CR- PS Trial Liner – #2 / 12mm          | 1    |
| L1   | 9066.35.214 | CR- PS Trial Liner – #2 / 14mm          | 1    |
| L1   | 9066.35.217 | CR- PS Trial Liner – #2 / 17mm          | 1    |
| L1   | 9066.35.220 | CR- PS Trial Liner – #2 / 20mm          | 1    |
| L1   | 9066.35.310 | CR- PS Trial Liner – #3 / 10mm          | 1    |
| L1   | 9066.35.312 | CR- PS Trial Liner – #3 / 12mm          | 1    |
| L1   | 9066.35.314 | CR- PS Trial Liner – #3 / 14mm          | 1    |
| L1   | 9066.35.317 | CR- PS Trial Liner – #3 / 17mm          | 1    |
| L1   | 9066.35.320 | CR- PS Trial Liner – #3 / 20mm          | 1    |
| L1   | 9066.35.410 | CR- PS Trial Liner – #4 / 10mm          | 1    |
| L1   | 9066.35.412 | CR- PS Trial Liner – #4 / 12mm          | 1    |
| L1   | 9066.35.414 | CR- PS Trial Liner – #4 / 14mm          | 1    |
| L1   | 9066.35.417 | CR- PS Trial Liner – #4 / 17mm          | 1    |
| L1   | 9066.35.420 | CR- PS Trial Liner – #4 / 20mm          | 1    |
| L1   | 9066.35.510 | CR- PS Trial Liner – #5 / 10mm          | 1    |
| L1   | 9066.35.512 | CR- PS Trial Liner – #5 / 12mm          | 1    |
| L1   | 9066.35.514 | CR- PS Trial Liner – #5 / 14mm          | 1    |
| L1   | 9066.35.517 | CR- PS Trial Liner – #5 / 17mm          | 1    |
| L1   | 9066.35.520 | CR- PS Trial Liner – #5 / 20mm          | 1    |
| N1   | 9066.35.610 | Plier for Trial Tibial Liner extractor  | 1    |
| O1   | 9069.10.275 | Pin for Trial Tibial Plates             | 4    |
|      | 9066.55.950 | Sterilizable Box                        | 1    |


# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Instrument set

### ▼ 9066.56.000 CCK-H Base Set n. 2



Upper tray



Lower tray


| Ref. | CODE        | DESCRIPTION                              | Qty. |
|------|-------------|------------------------------------------|------|
| A2   | 9066.12.010 | Sickle                                   | 1    |
| B2   | 9066.12.030 | Starting Reamer                          | 1    |
| C2   | 9066.15.050 | Distal Resection Guide                   | 1    |
| D2   | 9066.15.090 | Alignment Rod                            | 1    |
| E2   | 9066.22.120 | Tibial Cutting Guide                     | 1    |
| F2   | 9066.25.100 | Multifunction Impactor-Extractor         | 1    |
| G2   | 9066.25.110 | Tibial Impactor                          | 1    |
| H2   | 9066.25.160 | 0mm Tibial Depth Resection Gauge         | 1    |
| I2   | 9066.25.190 | Multifunction Extractor                  | 1    |
| J2   | 9066.30.020 | Awl                                      | 1    |
| K2   | 9066.30.100 | Tibial Winged Broach                     | 1    |
| L2   | 9066.30.110 | Guide for Winged Broach/ Reamer          | 1    |
| M2   | 9066.30.120 | Femoral Extractor/Positioner             | 1    |
| N2   | 9066.35.620 | Fiches Extractor Plier                   | 1    |
| O2   | 9066.35.621 | Inertial Beater                          | 1    |
| P2   | 9066.47.055 | Augment Tibial Resection Guide           | 1    |
| Q2   | 9066.47.064 | Joint line guide for femoral distal mask | 1    |
| R2   | 9066.47.065 | Joint line guide for femoral box mask    | 1    |
| S2   | 9066.47.225 | Tibial Reamer                            | 1    |
| T2   | 9066.47.650 | Trial Distal Augment - H 5mm             | 2    |
| T2   | 9066.47.655 | Trial Distal Augment - H 10mm            | 2    |
| U2   | 9066.50.130 | Osteotome #0                             | 1    |
| V2   | 9069.10.220 | Femoral Impactor                         | 1    |
| W2   | 9069.50.040 | Tibial Liner Extractor                   | 1    |
| X2   | 9095.10.337 | Flat Rasp                                | 1    |
|      | 9066.56.950 | Sterilizable Box                         | 1    |

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### Instrument set

#### ▼ 9066.41.000 Patellar Set n. 7


| Ref. | CODE        | DESCRIPTION                            | Qty. |
|------|-------------|----------------------------------------|------|
| A7   | 9066.40.050 | Patellar Gauge                         | 1    |
| B7   | 9066.40.010 | Patellar Pliers                        | 1    |
| C7   | 9066.40.020 | Patellar Vice                          | 1    |
| D7   | 9066.40.025 | Patellar Drilling Guide                | 1    |
| E7   | 9066.40.070 | Patellar Mask                          | 1    |
| F7   | 9066.40.060 | Patellar Drill                         | 1    |
| G7   | 9066.40.128 | Trial Patellar Dome - Dia. 28mm h 8mm  | 1    |
| G7   | 9066.40.132 | Trial Patellar Dome - Dia. 32mm h 8mm  | 1    |
| G7   | 9066.40.032 | Trial Patellar Dome - Dia. 32mm h 10mm | 1    |
| G7   | 9066.40.035 | Trial Patellar Dome - Dia. 35mm h 10mm | 1    |
| G7   | 9066.40.038 | Trial Patellar Dome - Dia. 38mm h 10mm | 1    |
| G7   | 9066.40.041 | Trial Patellar Dome - Dia. 41mm h 10mm | 1    |
|      | 9066.41.950 | Sterilizable Box                       | 1    |



# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Instrument Set

▼ 9066.51.000 CCK-H System Additional Set - Multigen Plus



| Ref. | CODE        | DESCRIPTION                        | Qty. |
|------|-------------|------------------------------------|------|
| A51  | 9066.15.105 | Pin Dia. 3.5x105mm                 | 3    |
| B51  | 9066.48.055 | IM Tibial Cutting Guide            | 1    |
| C51  | 9066.48.200 | Femoral Counter Torque             | 1    |
| D51  | 9066.48.205 | Tibial Pin Positioner              | 1    |
| E51  | 9066.48.210 | Tibial Pin Positioner Adapter      | 1    |
| F51  | 9066.48.220 | Femoral Box Gauge                  | 1    |
| G51  | 9095.10.226 | Dynamometric Screwdriver 3.5mm 5Nm | 1    |
|      | 9066.51.950 | Instrument Tray                    | 1    |

▼ 9095.11.751 Torque wrench

| Ref. | CODE        | DESCRIPTION   | Qty. |
|------|-------------|---------------|------|
|      | 9095.11.751 | Torque wrench | 1    |



**NOTE.** Space is forseen in the Instrument Tray 9066.51.950.

## REVISION KNEE SYSTEM SURGICAL TECHNIQUE

### Instruments needed for Hinge only


▼ Instruments for Hinge only (part of the basic CCK sets)

| Ref. | CODE        | DESCRIPTION             | Qty. |
|------|-------------|-------------------------|------|
| F47  | 9066.47.095 | H Winged Broach Guide   | 1    |
| E9   | 9066.48.050 | Im Cutting Guide H      | 1    |
| J9   | 9066.48.110 | H Trial Tibial Plate #1 | 1    |
| J9   | 9066.48.120 | H Trial Tibial Plate #2 | 1    |
| J9   | 9066.48.130 | H Trial Tibial Plate #3 | 1    |
| J9   | 9066.48.140 | H Trial Tibial Plate #4 | 1    |
| J9   | 9066.48.150 | H Trial Tibial Plate #5 | 1    |
| K9   | 9066.48.180 | H Tibial Keel           | 1    |

**NOTE.** *These instruments, together with the set 9066.51.000 CCK- H System Additional set and the torque wrench 9095.11.751 , are not needed for CCK procedures.*


# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Product combinations



# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Product codes



### ▼ CCK - CEMENTED FEMORAL COMPONENTS

CoCrMo

| Size | REF         |
|------|-------------|
| #1   | 6615.09.010 |
| #2   | 6615.09.020 |
| #3   | 6615.09.030 |
| #4   | 6615.09.040 |
| #5   | 6615.09.050 |



### ▼ H - CEMENTED FEMORAL COMPONENTS

CoCrMo+ CFR Peek

| Size | REF         |
|------|-------------|
| #1   | 6616.09.010 |
| #2   | 6616.09.020 |
| #3   | 6616.09.030 |
| #4   | 6616.09.040 |
| #5   | 6616.09.050 |



### ▼ CEMENTED FIXED TIBIAL PLATES

Ti6Al4V

| Size | REF         |
|------|-------------|
| #1   | 6624.15.110 |
| #2   | 6624.15.120 |
| #3   | 6624.15.130 |
| #4   | 6624.15.140 |
| #5   | 6624.15.150 |



### ▼ H - CEMENTED TIBIAL PLATES

CoCrMo+ CFR Peek

| Size | REF         |
|------|-------------|
| #1   | 6620.09.110 |
| #2   | 6620.09.120 |
| #3   | 6620.09.130 |
| #4   | 6620.09.140 |
| #5   | 6620.09.150 |



### ▼ CCK - TIBIAL LINERS WITH LOCKING SCREW

UHMWPE + Ti6Al4V

#### FOR FIXED TIBIAL PLATE #1

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #1   | 6645.50.110 | 10 mm     |
| #1   | 6645.50.112 | 12 mm     |
| #1   | 6645.50.114 | 14 mm     |
| #1   | 6645.50.117 | 17 mm     |
| #1   | 6645.50.120 | 20 mm     |

#### FOR FIXED TIBIAL PLATE #2

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #2   | 6645.50.210 | 10 mm     |
| #2   | 6645.50.212 | 12 mm     |
| #2   | 6645.50.214 | 14 mm     |
| #2   | 6645.50.217 | 17 mm     |
| #2   | 6645.50.220 | 20 mm     |

#### FOR FIXED TIBIAL PLATE #3

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #3   | 6645.50.310 | 10 mm     |
| #3   | 6645.50.312 | 12 mm     |
| #3   | 6645.50.314 | 14 mm     |
| #3   | 6645.50.317 | 17 mm     |
| #3   | 6645.50.320 | 20 mm     |

#### FOR FIXED TIBIAL PLATE #4

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #4   | 6645.50.410 | 10 mm     |
| #4   | 6645.50.412 | 12 mm     |
| #4   | 6645.50.414 | 14 mm     |
| #4   | 6645.50.417 | 17 mm     |
| #4   | 6645.50.420 | 20 mm     |

#### FOR FIXED TIBIAL PLATE #5

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #5   | 6645.50.510 | 10 mm     |
| #5   | 6645.50.512 | 12 mm     |
| #5   | 6645.50.514 | 14 mm     |
| #5   | 6645.50.517 | 17 mm     |
| #5   | 6645.50.520 | 20 mm     |

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Product codes



### ▼ H - TIBIAL LINERS WITH HINGE

UHMWPE + CoCrMo + Ti6Al4V

#### FOR H FEMORAL COMPONENT #1

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #1   | 6660.50.110 | 10 mm     |
| #1   | 6660.50.112 | 12 mm     |
| #1   | 6660.50.114 | 14 mm     |
| #1   | 6660.50.117 | 17 mm     |
| #1   | 6660.50.120 | 20 mm     |
| #1   | 6660.50.124 | 24 mm     |

#### FOR H FEMORAL COMPONENT #2

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #2   | 6660.50.210 | 10 mm     |
| #2   | 6660.50.212 | 12 mm     |
| #2   | 6660.50.214 | 14 mm     |
| #2   | 6660.50.217 | 17 mm     |
| #2   | 6660.50.220 | 20 mm     |
| #2   | 6660.50.224 | 24 mm     |

#### FOR H FEMORAL COMPONENT #3

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #3   | 6660.50.310 | 10 mm     |
| #3   | 6660.50.312 | 12 mm     |
| #3   | 6660.50.314 | 14 mm     |
| #3   | 6660.50.317 | 17 mm     |
| #3   | 6660.50.320 | 20 mm     |
| #3   | 6660.50.324 | 24 mm     |

#### FOR H FEMORAL COMPONENT #4

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #4   | 6660.50.410 | 10 mm     |
| #4   | 6660.50.412 | 12 mm     |
| #4   | 6660.50.414 | 14 mm     |
| #4   | 6660.50.417 | 17 mm     |
| #4   | 6660.50.420 | 20 mm     |
| #4   | 6660.50.424 | 24 mm     |

#### FOR H FEMORAL COMPONENT #5

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #5   | 6660.50.510 | 10 mm     |
| #5   | 6660.50.512 | 12 mm     |
| #5   | 6660.50.514 | 14 mm     |
| #5   | 6660.50.517 | 17 mm     |
| #5   | 6660.50.520 | 20 mm     |
| #5   | 6660.50.524 | 24 mm     |



### ▼ CCK/H - FEMORAL MODULES Ti6Al4V

| REF         | LENGTH | OFFSET          |
|-------------|--------|-----------------|
| 6647.15.410 | SHORT  | R-L Straight    |
| 6647.15.420 | SHORT  | R +3mm - L -3mm |
| 6647.15.430 | SHORT  | R -3mm - L +3mm |
| 6647.15.440 | LONG   | R-L Straight    |
| 6647.15.450 | LONG   | R +3mm - L -3mm |
| 6647.15.460 | LONG   | R -3mm - L +3mm |



### ▼ CCK - TIBIAL MODULE Ti6Al4V

| REF         | LENGTH | OFFSET   |
|-------------|--------|----------|
| 6647.15.310 | SHORT  | Straight |
| 6647.15.320 | SHORT  | +3 mm    |
| 6647.15.330 | SHORT  | +6 mm    |
| 6647.15.340 | LONG   | Straight |
| 6647.15.350 | LONG   | +3 mm    |
| 6647.15.360 | LONG   | +6 mm    |



### ▼ H - TIBIAL MODULE Ti6Al4V

| REF         | LENGTH | OFFSET   |
|-------------|--------|----------|
| 6648.15.310 | SHORT  | Straight |
| 6648.15.320 | SHORT  | +3 mm    |
| 6648.15.330 | SHORT  | +6 mm    |
| 6648.15.340 | LONG   | Straight |
| 6648.15.350 | LONG   | +3 mm    |
| 6648.15.360 | LONG   | +6 mm    |



### ▼ CCK/H - STEM Ti6Al4V

| REF         | DIAMETER | LENGTH |
|-------------|----------|--------|
| 6647.15.140 | 14 mm    | 60 mm  |
| 6647.15.160 | 16 mm    | 60 mm  |
| 6647.15.180 | 18 mm    | 85 mm  |
| 6647.15.200 | 20 mm    | 85 mm  |
| 6647.15.220 | 22 mm    | 110 mm |
| 6647.15.240 | 24 mm    | 110 mm |

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Product codes



### ▼ PATELLAR PROSTHESES

#### UHMWPE

|   | REF         | DIAMETER | THICKNESS |
|---|-------------|----------|-----------|
| ■ | 6695.50.103 | 28 mm    | 8 mm      |
| ■ | 6695.50.105 | 32 mm    | 8 mm      |
| ■ | 6695.50.005 | 32 mm    | 10 mm     |
| ■ | 6695.50.010 | 35 mm    | 10 mm     |
| ■ | 6695.50.020 | 38 mm    | 10 mm     |
| ■ | 6695.50.030 | 41 mm    | 10 mm     |

■ upon request



### ▼ TIBIAL AUGMENTS

#### Ti6Al4V

| FOR FIXED TIBIAL PLATE #1 |             |           |
|---------------------------|-------------|-----------|
| Size                      | REF         | THICKNESS |
| #1                        | 6625.15.010 | 7 mm      |
| #1                        | 6625.15.020 | 12 mm     |
| FOR FIXED TIBIAL PLATE #2 |             |           |
| Size                      | REF         | THICKNESS |
| #2                        | 6626.15.010 | 7 mm      |
| #2                        | 6626.15.020 | 12 mm     |
| FOR FIXED TIBIAL PLATE #3 |             |           |
| Size                      | REF         | THICKNESS |
| #3                        | 6628.15.010 | 7 mm      |
| #3                        | 6628.15.020 | 12 mm     |
| FOR FIXED TIBIAL PLATE #4 |             |           |
| Size                      | REF         | THICKNESS |
| #4                        | 6630.15.010 | 7 mm      |
| #4                        | 6630.15.020 | 12 mm     |
| FOR FIXED TIBIAL PLATE #5 |             |           |
| Size                      | REF         | THICKNESS |
| #5                        | 6632.15.010 | 7 mm      |
| #5                        | 6632.15.020 | 12 mm     |



## ▼ H - TIBIAL AUGMENTS

Ti6Al4V

### FOR H TIBIAL PLATE #1

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #1   | 6621.15.107 | 7 mm      |
| #1   | 6621.15.112 | 12 mm     |

### FOR H TIBIAL PLATE #2

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #2   | 6621.15.207 | 7 mm      |
| #2   | 6621.15.212 | 12 mm     |

### FOR H TIBIAL PLATE #3

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #3   | 6621.15.307 | 7 mm      |
| #3   | 6621.15.312 | 12 mm     |

### FOR H TIBIAL PLATE #4

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #4   | 6621.15.407 | 7 mm      |
| #4   | 6621.15.412 | 12 mm     |

### FOR H TIBIAL PLATE #5

| Size | REF         | THICKNESS |
|------|-------------|-----------|
| #5   | 6621.15.507 | 7 mm      |
| #5   | 6621.15.512 | 12 mm     |

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE


## Product codes



### ▼ DISTAL FEMORAL AUGMENTS

Ti6Al4V

| FOR FEMORAL COMPONENT #1 |             |           |
|--------------------------|-------------|-----------|
| Size                     | REF         | THICKNESS |
| #1                       | 6671.15.105 | 5 mm      |
| FOR FEMORAL COMPONENT #2 |             |           |
| Size                     | REF         | THICKNESS |
| #2                       | 6671.15.205 | 5 mm      |
| FOR FEMORAL COMPONENT #3 |             |           |
| Size                     | REF         | THICKNESS |
| #3                       | 6671.15.305 | 5 mm      |
| #3                       | 6671.15.310 | 10 mm     |
| FOR FEMORAL COMPONENT #4 |             |           |
| Size                     | REF         | THICKNESS |
| #4                       | 6671.15.405 | 5 mm      |
| #4                       | 6671.15.410 | 10 mm     |
| FOR FEMORAL COMPONENT #5 |             |           |
| Size                     | REF         | THICKNESS |
| #5                       | 6671.15.505 | 5 mm      |
| #5                       | 6671.15.510 | 10 mm     |



▼ POSTERIOR FEMORAL AUGMENTS

Ti6Al4V

| FOR FEMORAL COMPONENT #1 |             |           |
|--------------------------|-------------|-----------|
| Size                     | REF         | THICKNESS |
| #1                       | 6672.15.105 | 5 mm      |
| FOR FEMORAL COMPONENT #2 |             |           |
| Size                     | REF         | THICKNESS |
| #2                       | 6672.15.205 | 5 mm      |
| FOR FEMORAL COMPONENT #3 |             |           |
| Size                     | REF         | THICKNESS |
| #3                       | 6672.15.305 | 5 mm      |
| #3                       | 6672.15.310 | 10 mm     |
| FOR FEMORAL COMPONENT #4 |             |           |
| Size                     | REF         | THICKNESS |
| #4                       | 6672.15.405 | 5 mm      |
| #4                       | 6672.15.410 | 10 mm     |
| FOR FEMORAL COMPONENT #5 |             |           |
| Size                     | REF         | THICKNESS |
| #5                       | 6672.15.505 | 5 mm      |
| #5                       | 6672.15.510 | 10 mm     |

# REVISION KNEE SYSTEM SURGICAL TECHNIQUE

## Tables summary

Here following are listed some useful tables to be checked before opening the definitive implants, such as the femoral or tibial modules:

### ▼ FEMORAL OR TIBIAL MODULE LENGHT

| REAMER REFERENCE   | TRIAL MODULE           | DEFINITIVE MODULE |
|--------------------|------------------------|-------------------|
| <b>S Reference</b> | S Module               | S Module          |
| <b>L Reference</b> | S Module + L Extension | L Module          |

### ▼ TIBIAL OFFSET

| COMPASS   | TIBIAL TRIAL MODULE | TIBIAL DEFINITIVE MODULE |
|-----------|---------------------|--------------------------|
| <b>0</b>  | Straight Module     | Straight Module          |
| <b>+3</b> | +3 Module           | +3 Module                |
| <b>+6</b> | +6 Module           | +6 Module                |

### ▼ FEMORAL OFFSET

| ALIGNER           | FEMORAL TRIAL MODULE | FEMORAL DEFINITIVE MODULE |
|-------------------|----------------------|---------------------------|
| <b>R-L</b>        | R/L Module           | R/L Module                |
| <b>R+3 or L-3</b> | R+3/L-3 Module       | R+3/L-3 Module            |
| <b>R-3 or L+3</b> | R-3/L+3 Module       | R-3/L+3 Module            |

Here following are listed the total lengths of the CCK/H stems assembled with femoral or with tibial modules:

### ▼ CCK/H STEMS AND FEMORAL MODULES ASSEMBLY LENGTH

| CCK/H PRESS FIT<br>STEM DIAMETERS | F1                              |                                |
|-----------------------------------|---------------------------------|--------------------------------|
|                                   | STEM +<br>SHORT FEMORAL MODULE* | STEM +<br>LONG FEMORAL MODULE* |
| Dia. 14 mm                        | 91 mm                           | 116 mm                         |
| Dia. 16 mm                        | 91 mm                           | 116 mm                         |
| Dia. 18 mm                        | 116 mm                          | 141 mm                         |
| Dia. 20 mm                        | 116 mm                          | 141 mm                         |
| Dia. 22 mm                        | 141 mm                          | 166 mm                         |
| Dia. 24 mm                        | 141 mm                          | 166 mm                         |


\* for the femoral component 18,5 mm must be added to obtain the total length from the femoral distal resection to the end of the stem.

### ▼ CCK/H STEMS AND TIBIAL MODULES ASSEMBLY LENGTH

| CCK/H PRESS FIT<br>STEM DIAMETERS | T1                              |                                |
|-----------------------------------|---------------------------------|--------------------------------|
|                                   | STEM +<br>SHORT TIBIAL MODULE** | STEM +<br>LONG TIBIAL MODULE** |
| Dia. 14 mm                        | 87 mm                           | 112 mm                         |
| Dia. 16 mm                        | 87 mm                           | 112 mm                         |
| Dia. 18 mm                        | 112 mm                          | 137 mm                         |
| Dia. 20 mm                        | 112 mm                          | 137 mm                         |
| Dia. 22 mm                        | 137 mm                          | 162 mm                         |
| Dia. 24 mm                        | 137 mm                          | 162 mm                         |

\*\* for the tibial component 22 mm must be added to obtain the total length from the tibial resection to the end of the stem.

\*\*\* for the individual lengths of the modules, please refer to page 87.







**Limacorporate S.p.A.**

Via Nazionale, 52  
33038 Villanova di San Daniele del Friuli  
Udine - Italy  
T +39 0432 945511  
F +39 0432 945512  
info@limacorporate.com  
limacorporate.com

**Lima Implantas sl**

Calle Asura n. 97  
Madrid 28043  
España

**Lima France sas**

1, Allée des Alisiers  
Immeuble le Galilée  
69500 Bron  
France  
T +33 4 87 25 84 30  
F +33 4 42 04 17 25  
info@limafrance.com

**Lima O.I. doo**

Ante Kovacica, 3  
10000 Zagreb - Croatia  
T +385 (0) 1 2361 740  
F +385 (0) 1 2361 745  
lima-oi@lima-oi.hr

**Lima Switzerland sa**

Birkenstrasse, 49  
CH-6343 Rotkreuz - Zug  
Switzerland  
T +41 (0) 41 747 06 60  
F +41 (0) 41 747 06 69  
info@lima-switzerland.ch

**Lima Japan kk**

Shinjuku Center Building, 29<sup>th</sup> floor  
1-25-1, Nishi-shinjuku, Shinjuku,  
Tokyo 163-0629 - Japan  
T +81 3 5322 1115  
F +81 3 5322 1175

**Lima CZ sro**

Do Zahrádek I., 157/5  
155 21 Praha 5 - Zličín  
Czech Republic  
T +420 222 720 011  
F +420 222 723 568  
info@limacz.cz

**Lima Deutschland GmbH**

Kapstadttring 10  
22297 Hamburg - Germany  
T +49 40 6378 4640  
F +49 40 6378 4649  
info@lima-deutschland.com

**Lima Austria GmbH**

Seestadtstrasse 27 / Top 6-7  
1220 Wien - Austria  
T +43 (1) 2712469  
F +43 (1) 2712469101  
office@lima-austria.at

**Lima SK s.r.o.**

Cesta na štadión 7  
974 04 Banská Bystrica - Slovakia  
T +421 484 161 126  
F +421 484 161 138  
info@lima-sk.sk

**Lima Netherlands**

Havenstraat 30  
3115 HD Schiedam  
The Netherlands  
T +31 (0) 10 246 26 60  
F +31 (0) 10 246 26 61  
info@limanederland.nl  
limanederland.nl

**Lima Implantas Portugal S.U. Lda**

Rua Olavo D'Eça Leal N°6 Loja-1  
1600-306 Lisboa - Portugal  
T +35 121 727 233 7  
F +35 121 296 119 2  
lima@limaportugal.com

**Lima Orthopaedics Australia Pty Ltd**

Unit 1, 40 Ricketts Rd  
Mt Waverley 3149  
Victoria Australia  
T +61 (03) 9550 0200  
F +61 (03) 9543 4003  
limaortho.com.au

**Lima Orthopaedics New Zealand Ltd**

20 Crummer Road  
Auckland 1021  
New Zealand  
T +64 93606010  
F +64 93606080

**Lima Orthopaedics UK Limited**

Unit 1, Campus 5  
Third Avenue  
Letchworth Garden City  
Herts, SG6 2JF  
United Kingdom  
T +44 (0) 844 332 0661  
F +44 (0) 844 332 0662

**Lima USA Inc.**

2001 NE Green Oaks Blvd., Suite 100  
Arlington, TX 76006  
T +1 817-385-0777  
F +1 817-385-0377

**Lima Sweden AB**

Företagsallén 14 B  
SE-184 40 ÅKERSBERGA  
Sweden  
T +46 8 544 103 80  
F +46 8 540 862 68  
www.links sweden.se

**Lima Italy**

Centro Direzionale Milanofiori  
Strada 1 - Palazzo F9  
20057 Assago - Milano - Italy  
T +39 02 57791301

**Lima Korea Co. Ltd**

11 FL., Zero Bldg.  
14 Teheran Road 84 GLL  
Gangnam Gu, Seoul 135-845, South Korea  
T +82 2 538 4212  
F +82 2 538 0706

**Lima do Brasil EIRELI**

Al. Campinas, 728, second floor,  
rooms 201, 202, 203 and 204,  
Edifício Engenheiro Antonio Silva,  
Zip Code 01404-001, in the City of São Paulo,  
State of São Paulo  
Brasil

**Lima Belgium sprl**

Chaussée de Wavre 504, bte 48  
1390 Grez-Doiceau - Belgium  
T +32 (0) 10 888 804  
F +32 (0) 10 868 117  
info@limabelgium.be

**Lima Denmark ApS**

Lyngbækgårds Allé 2  
2990 Nivå - Denmark  
T +45 45860028  
F +45 4586 0068  
mail@Lima-Denmark.dk

**Lima Turkey Ortopedi A.S.**

Serifali Mah. Hendem CD. Canan  
Residence No: 54/C D:2 OFIS-A2,  
34775 Umraniye / Istanbul  
Turkey  
T +90 (216) 693 1373  
F +90 (216) 693 2212  
info@lima-turkey.com.tr

**Lima Orthopaedics South Africa**

Northlands Deco Park, Stand 326  
10 New Market street  
Design Boulevard  
Northriding  
2189

**Lima Polska Sp. z o.o.**

Ul. Łopuszańska 95  
02-457 Warszawa  
Poland  
T 0048 22 6312786  
F 0048 22 6312604  
biuro@limapolska.pl

**励玛（北京）医疗器械有限公司**

**Lima (Beijing) Medical Devices Co., Ltd.**

中国北京市朝阳区利泽中二路1号中辰大厦6层 616室  
Room 616, 6/F Zhongchen Building, No.1 Lize Zhong 2 Road  
Chaoyang District,  
Beijing - PR China  
limachinaoffice@limacorporate.com

The people appearing in the photographs on this publication are models and used for illustrative purposes only.

Prescription use only: sold only on or by the order of a physician.

This publication is not intended for distribution in the U.S.

B.6615.23.051.1

052100



limacorporate.com

