

Chirurgia del piede / *Foot surgery*

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Foot surgery

Chirurgia del piede

Product description
& surgical technique

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
Clinical cases

Totally conservative surgery on advanced rheumatoid forefoot:

- osteotomy of the 1st phalanx of the hallux by shape memory staple
- osteotomy of the 1st metatarsus by two double threaded screws
- osteotomy of the lateral metatarsus by self-breaking screws.

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
Indications

Ø 3 mm HVS SCREW:

- first metatarsal correction
- fixation of intra-articular and peri-articular fragments.

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Materials

Titanium alloy Ti6Al4V in conformity to ISO 5832-3.

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Technical features

- Self-drilling and self-tapping.
- Square coupling socket with screwdriver.
- Double threaded titanium cannulated screws with differentiated pitch are suitable in case of fractures where the lack of tissue makes the use of the traditional cancellous screw inappropriate or where a fragments alignment is needed.
- The headless screw must be completely embedded in the bone thus avoiding any protrusion which could cause tissue irritations.
- Its special thread assures fracture compression and stability.

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**HVS self-drilling and self-tapping
 double threaded cannulated screw**
**Vite a doppia filettatura cannulata
 autoperforante e autofilettante HVS**

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


Surgical technique


1. Lateral release and exostosis resection.

2. Scarf osteotomy.

3. Kirschner wire insertion and measurement with HME 353-01 measuring device.

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


Indications

SELF-BREAKING SCREW:

- diaphyseal fractures of phalanges and joints
- fractures with small-sized bone fragments
- osteotomy of lateral metatarsus.

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

Technical features

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The self-breaking screws are suitable for hand and foot surgery.

The diameter of these screws is 2 mm or 2.7 mm while the cylindrical screw holder stem, to be inserted in a Jacob's chuck, has a $\varnothing$ 1.6 mm.

Peculiarities:

- fast insertion through a powered drill
- self-drilling and self-tapping
- wide range of sizes
- type 3 anodic oxidation for color code differentiation of diameters.

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
Product codes

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Self-breaking screws Viti autorompenti		
Ti6Al4V		
Length Lunghezza mm	Ø 2 mm	Ø 2.7 mm
7	HMV 390-07-1ST	-
9	HMV 390-09-1ST	-
10	HMV 390-10-1ST	-
11	HMV 390-11-1ST	HMV 395-11-1ST
12	HMV 390-12-1ST	HMV 395-12-1ST
13	HMV 390-13-1ST	HMV 395-13-1ST
14	HMV 390-14-1ST	HMV 395-14-1ST
15	HMV 390-15-1ST	HMV 395-15-1ST
16	-	HMV 395-16-1ST
17	HMV 390-17-1ST	HMV 395-17-1ST
18	-	HMV 395-18-1ST
19	HMV 390-19-1ST	HMV 395-19-1ST
20	-	HMV 395-20-1ST
22	-	HMV 395-22-1ST

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

Surgical technique

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
1. Prepare a pre-hole with a $\varnothing$ 1 mm Kirschner wire if necessary, then place the screw on the external cortex and turn on the drill by pushing to perforate the cortex.

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
2. Avoid sudden movements while inserting the screw to prevent the early breaking of the head. In this case a special screwdriver (HMTV 392) is available to complete the screwing.

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3. The detachment of the barrel occurs during compression as soon as the screw head reaches the cortex level.

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4. Remove the excess bone section to restore the good joint.

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Indications

SERRATED STAPLE:

- Scaphoid fractures
- Hand and foot surgery

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Materials

The device is manufactured in AISI 316L ESR stainless steel in conformity to ISO 5832-1.

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Product codes

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Serrated straight staple Cambra dritta dentata

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


Surgical technique

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
1. Prepare the site of the compressor bits by making holes with a $\varnothing$ 1 mm wire or bit; apply the compressor HMA 246-1 once osteotomy has been made.

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
2. Apply the mask HMA 246-4 to evaluate the exact size of the staple to be used, drill through the mask holes by using a $\varnothing$ 1 mm wire or bit to create the site for the staple legs.

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3. Staple insertion with the proper staple holder impactor HMA 246-2 or HMA 246-6.

3. Applicare la cambra con l'apposito impattatore porta-cambra HMA 246-2 o HMA 246-6.

4. Use the final impactor HMA 246-3 or HMA 246-5 to definitively hit the staple. Remove the compressor HMA 246-1.

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Indications

HOT MEMORY STAPLES

Osteosynthesis for treatment of small bone osteotomies, arthrodesis and internal fixation of fractures.

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Nitinol conforme alla norma ASTM F2063.

Technical features

The hot memory staples are made of a special nickel-titanium alloy with thermoelastic transformation.

The use of electrosurgical knife is recommended in clot position. Power intensity adjustment allows the device to regain the original shape (memory effect) set between temperature range of 43° - 55°C. It is above 43°C that heating flow activates the material transition.

The staple legs closure and therefore the compression is adjustable with the heat.

Staples are supplied in sterile packaging in a suitable support that avoids the deformation in case of accidental heat sources. They have to be kept at room temperature.

The staple section is either square (e.g. 1.2x1.2 mm, 1.5x1.5 mm) or rectangular (e.g. 3x2 mm).

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Surgical technique

The holes must have a larger diameter than the staple one (e.g. $\varnothing$ 2 mm for staple with section 1.5 mm, $\varnothing$ 1.5 mm for staple with section 1.2 mm and $\varnothing$ 3.5 mm for staple with section 3 mm).

It is advisable to adjust the hole guide 2 mm longer than bridge width (e.g. staple HMA B 12-13-15 has to be adjusted at 14 mm) in order to obtain a better compression.

Use the forceps to ease the insertion of the two legs and then use a straight hammer.

TO OBTAIN BETTER COMPRESSION WITH ASYMMETRIC STAPLE

1. Set the drill guide 2 mm larger than the staple's width (eg 12 mm than set 14 mm)
2. Perform drilling on both holes
3. Insert the longer leg in the proximal hole and then insert the shorter into the distal hole by 'forcing' the staple's leg inside (use a klemmer)
4. Push the staple's bridge close to the medial cortex
5. Heat for 3/5 seconds with a soft sliding up and down on the staples's bridge

WARNINGS

USE THE ELECTROSURGICAL KNIFE FOR THE NECESSARY STAPLE DEFORMATION TIME ONLY. A PROLONGED USE COULD DAMAGE THE SURROUNDING TISSUES.

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
Surgical technique

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
OSTEOTOMY OF THE FIRST HALLUX PHALANX

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
1. Use a $\varnothing$ 2 mm drill bit at 8-10 mm from the spacing in horizontal direction and parallel to the articular surface to perform a proximal pilot hole. A wire of the same diameter should be inserted in the hole as a reference point. Measure the pilot hole depth in order to establish the length of the longer leg.

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


Surgical technique

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5. Distal hole preparation. Once the guide has been set on the size of the staple to be used, position the guide on the proximal wire which is already inserted and make the second distal hole with a $\varnothing$ 2 mm drill bit therefore.

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
6. Position the staple by using the specific staple holder and use the electrosurgical knife in order to allow the device to regain the original shape (memory effect).

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VITE CALCANEUS STOP

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Materials

Ti6Al4V titanium alloy in conformity to ISO 5832-3.

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Technical features

This calcaneal endorthesis has been designed to reach an extra-articular arthrorisis of coxa pedis, able to compensate ratios in case of astragalic protrusion and normalize the plantar helix in case of levelling dysfunction. The goal is the stable reduction of the protusion by surgery. The screw is self-tapping and has a tapered head eccentrically positioned towards the thread axis. The technique is based on the following principles:

- repositioning of the astragalus on the calcaneus by abducting the first and varising the second one, in order to adjust the ratio and restore the physiological evolutive lines;
- artificial maintenance of the new articular position up to the self-support allowed by the augmentative remodelling.

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Product codes

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ø 7 mm cannulated calcaneo stop screw Vite calcaneo stop cannulata ø 7 mm

Ti6Al4V

Length Lunghezza mm	
28	HMV 445-28ST
32	HMV 445-32ST
36	HMV 445-36ST

Sterile single packaging.

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Instrument set

Strumentario


Instrument set Strumentario


Ref.	Description	Descrizione	
1	ø 2 mm guide wire length 230 mm	Filo guida ø 2 mm lunghezza 230 mm	HMV 353-04
2	ø 4.5 mm cannulated drill bit	Punta cannulata ø 4,5 mm	HMV 650-000
3	Cannulated screwdriver hex 3.5	Cacciavite cannulato esagono 3,5	HMV 941
4	Measuring device	Misuratore	HM 00110
-	Tray for instrument	Vassoio per strumenti	HM 106-000

Surgical technique

The technique basically consists of the manual or instrumental reduction of the inter-astragal-calcaneus detorsion and the maintenance of the adequate ratio through the screw.

1. Insertion of a $\varnothing$ 2 mm guide wire and evaluation of the screw size to be used.

2. Use the $\varnothing$ 4.5 mm cannulated drill bit to perforate the tarsal sinus floor (optional step).

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
Surgical technique

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3. Implant the screw in the anterior process of the calcaneus, with dorso-plantar, lateral-medial and posterior-anterior direction, up to its anterior positioning to the lateral astragalic process. The self-tapping and slightly tapered tip allows an easier insertion.

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At the end of surgery eccentricity of the tapered head towards the screw axis allows an easier procedure to bring easily the calcaneous axis on the leg one, by screwing clockwise and counterclockwise.

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Herbert TJ, Fisher WE.
Management of the fractured scaphoid using a
new bone screw
J Bone Joint Surg Br. 1984 Jan;66(1):114-23.

Nery C, Barroco R, Réssio C.
Biplanar chevron osteotomy
Foot Ankle Int. 2002 Sep;23(9):792-8.

Ramaswamy R, Evans S, Kosashvili Y.
Holding power of variable pitch screws in osteoporotic,
osteopenic and normal bone: Are all screws created equal?
Injury. 2010 Feb;41(2):179-83.

Self-breaking screws Viti autorompenti

Reinhard S, Trnka HJ.
Metatarsalgia: Distal Metatarsal Osteotomies
Foot Ankle Clin N Am 16 (2011) 583–595.

Trnka HJ, Gebhard C, Mühlbauer M, Ivanic G, Ritschl P.
The Weil osteotomy for treatment of dislocated lesser
metatarsophalangeal joints: good outcome in 21 patients
with 42 osteotomies
Acta Orthop Scand. 2002 Apr;73(2):190-4.

Barouk LS, Barouk P.
Joint-preserving surgery in rheumatoid forefoot: preliminary
study with more-than-two-year follow-up
Foot Ankle Clin. 2007 Sep;12(3):435-54, vi.

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Shabalovskaya SA.
Surface, corrosion and biocompatibility aspects of Nitinol as an
implant material
Biomed Mater Eng. 2002;12(1):69-109.

Mereau TM, Ford TC.
Nitinol compression staples for bone fixation in foot surgery
J Am Podiatr Med Assoc. 2006 Mar-Apr;96(2):102-6.

Bensmann G, Baumgart F, Haaster J.
Nickel-titanium osteosynthesis clips. Advantages of new implant
confirmed in clinical application
Medical Focus, Exportmarkt-Organisation,
Vogel-Verlag Würzburg, Heft 3/1983.

Castleman LS, Motzkin SM, Alicandri FP, Bonawit VL., Johnson AA
Biocompatibility of nitinol alloy as an implant material
J Biomed Mater Res. 1976 Sep; 10(5): 695-731.

Calcaneo stop HMV 445

Kellermann P, Roth S, Gion K, Boda K, Tóth K.
Calcaneo-stop procedure for paediatric flexible flatfoot
Arch Orthop Trauma Surg. 2011 Oct;131(10):1363-7.

Richter M, Zech S.
Arthrosis with calcaneostop screw in children corrects
Talo-1st Metatarsal-
Index (TMT-Index)
Foot Ankle Surg. 2013 Jun;19(2):91-5.

Jerosch J, Schunck J, Abdel-Aziz H.
The stop screw technique-a simple and reliable method
in treating flexible flatfoot in children
Foot Ankle Surg. 2009;15(4):174-8.

[illegible]

Info

LSM-MED S.r.l. reserves
the right to make changes.

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
*Implants and instruments are manufactured by LSM-MED S.r.l.
and distributed and marketed by Limacorporate S.p.A.*

Distribuito da / *Distributed by:*

Limacorporate S.p.A.

Via Nazionale, 52
33038 Villanova di San Daniele del Friuli
Udine - Italy
T +39 0432 945511
F +39 0432 945512
info@limacorporate.com
limacorporate.com

B.0042.32.000.1

022000

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OVERMED S.r.l.
Impianti e strumenti sono prodotti da LSM-MED S.r.l.
e distribuiti e commercializzati da Limacorporate S.p.A.