



Highly
complex
glenoid



Highly
complex
elbow



Highly
complex
acetabulum



PROMADE

Bespoke Design Service

+ DESIGN SERVICE


*to develop patient specific implants
based on 3D printing know-how.*

+ BESPOKE SOLUTIONS

to manage the most challenging cases.

+ CREATING BETTER TOGETHER

*Modular design to guarantee
compatibility with full
Enovis portfolio.*



DESIGN SERVICE

- **Dedicated Engineer for complete case management and development.**
Direct partnership with your engineer to define design, surgical approach, surgical technique. Your engineer is accountable for the entire case from beginning to shipment.
- **Models and PSIs.** ProMade Service provides everything required to carefully plan the case outside the OR and seamlessly execute it during surgery.
- **Fully sterile Implant and PSIs** to reduce infection risks.




CONFIDENCE THROUGH PRIMARY STABILITY

- Increased **primary stability** thanks to **Trabecular Titanium** high friction between bone and implant to reduce micro-motion risk.^[1-3]
- **TT** low elastic modulus allows proper **load transfer** to the bone and reduces stress shielding and bone resorption.^[1-3]
- A defect filling approach to implant design, with irregular shaped bone interfaces, aims to provide a **highly stable initial placement**.


FULLY PLANNED BIOMECHANICAL RESTORATION

- **A comprehensive surgical plan**, using imaging data and clinical insight, to fully achieve desired patient biomechanics.
- **Design versatility** to accurately restore patient's **anatomy** and filling of **defects**.
- **Expanded modularity** facilitated by modular connections with Enovis portfolio providing a wide range of surgical options.



LEGEND

- COR: Centre Of Rotation
- ALP: Acromial Lateral Point
- CBP: Coracoid Base Point



Bibliography



B.0013.83.022.1

042500

Limacorporate S.p.A.
Via Nazionale, 52
33038, Villanova di San Daniele del Friuli
Italy

This publication is not intended for distribution in the U.S.
This document is not intended for circulation, publication or distribution to the public.
It is supplied for the exclusive use of Healthcare Professionals.